SilentCrusher M

**Betriebsanleitung
Instruction Manual
Mode d'Emploi
Instrucciones de
Empleo
Istruzioni per l'uso**

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EN 61010

Thank you for purchasing this equipment. You have acquired a product manufactured and tested by Heidolph Instruments in accordance with DIN EN ISO 61010. With this device you will be able to perform your work flawlessly and problem-free.

CONTENTS

CONTENTS	30
COMPONENTS SUPPLIED AND ACCESSORIES.....	32
GENERAL INFORMATION	33
SAFETY INSTRUCTIONS.....	33
ASSEMBLY	34
1. SET UP STAND (IF ORDERED AS AN OPTION).....	34
2. MOUNT DRIVE UNIT	34
3. MOUNT VESSEL WITH VESSEL RETAINER (IF ORDERED AS AN OPTION)	34
4. INSERT DISPERSION TOOL (STANDARD)	35
5. INSERT SEALED DISPERSION TOOL	36
6. POSITIONING THE DISPERSION TOOL IN THE FLUID	37
7. OPERATING UNDER PRESSURE OR IN A VACUUM WITH THE SEALED DISPERSION TOOL	37
8. CONNECTING THE POWER SUPPLY (UNIT) (S)	37
SERVICE AND OPERATION	38
9. START UP DEVICE	38
9.1 <i>Speed limitation of the drive</i>	39
10. TOOL REPLACEMENT	39
11. LENGTH OF OPERATION	40
CLEANING AND MAINTENANCE	41
12. DRIVE UNIT.....	41
13. DISPERSION TOOL.....	42
13.2. <i>Cleaning:</i>	42
13.3. <i>Disassembly, Standard Dispersion Tool:</i>	42
13.4. <i>Disassembly, Sealed Dispersion Tool:</i>	43
13.5. <i>Sterilization of the Dispersion Tool:</i>	45
13.6. <i>Assembly of the Dispersion Tool:</i>	45
SPARE AND WEAR PARTS.....	45
DISASSEMBLY, TRANSPORT AND STORAGE	49
DISPOSAL.....	49
TROUBLESHOOTING.....	49
TECHNICAL DATA	50
WARRANTY, LIABILITY AND COPYRIGHTS.....	56
QUESTIONS / REPAIRS	57

Important Notices

Connection Instructions / Power Connection

Caution, must be observed

Caution, Fire or Explosion Hazard

Repair / Maintenance Instructions

Instructions for personal protection

Caution: Risk of Combustion, Hot Surfaces

COMPONENTS SUPPLIED AND ACCESSORIES




Components supplied

Description	Quantity	Part Number 100 - 230 V / 50/60 Hz
SilentCrusher M eu (European plug) or	1	595-06000-00
SilentCrusher M us (US plug) or	1	595-06000-01
SilentCrusher M ch (CH plug)	1	595-06000-04
Power Supply (unit)	1	11-001-595-10
European Power Connection cable or	1	14-007-003-81
US Power Connection cable	1	14-007-003-89
Screwdriver	1	04-06-01-21-13
Operating Manual	1	01-005-004-46
O-ring 23 x 2mm (spare part)	3	23-08-06-03-17



Accessories (Optional)

Description	Part Number
Dispersion Tool 6 F Including 2 PTFE Bearings and spare part set	596-06010-00
Dispersion Tool 8 F Including 2 PTFE Bearings and spare part set	596-08010-00
Dispersion Tool 12 F Including 2 PTFE Bearings and spare part set	596-12010-00
Dispersion Tool 12 G Including 2 PTFE Bearings and spare part set	596-12020-00
Dispersion Tool 18 F Including 2 PTFE Bearings and spare part set	596-18010-00
Dispersion Tool 18 G Including 2 PTFE Bearings and spare part set	596-18020-00
Dispersion Tool 8 DF Including 2 PTFE Bearings and spare part set	596-08110-00
Dispersion Tool 12 DF Including 2 PTFE Bearings and spare part set	596-12110-00
Dispersion Tool 12 DG Including 2 PTFE Bearings and spare part set	596-12120-00
Dispersion Tool 18 DF Including 2 PTFE Bearings and spare part set	596-18110-00
Dispersion Tool 18 DG Including 2 PTFE Bearings and spare part set	596-18120-00
Dispersion Tool 22 F Including 2 PTFE Bearings and spare part set	596-22010-00
Dispersion Tool 22 G Including 2 PTFE Bearings and spare part set	596-22020-00
Plate Stand SilentCrusher M	570-17000-00
Vessel Retainer SilentCrusher M	570-30100-00
Sleeve coupling SilentCrusher M	570-31100-00

GENERAL INFORMATION

-  Please unpack the device carefully.
 Check for any damage and report damage or missing parts to the supplier immediately.
-  Read the operating instructions thoroughly and attentively and ensure that every operator of the device has read the operating instructions carefully before start-up.
-  Please keep the operating instructions in a location accessible to everyone.
-  The device is equipped with a EUROPEAN plug in accordance with the standards (DIN 49441 EEC 7/VII 10/ 16 A 250 V).
 For North America with a US STANDARD plug (NEMA Pub. No. WDI.1961 ASA C 73.1. 1961 page 8 15A 125V).
-  If you want to use the device in a country with a different outlet system, you must use an adapter or have the plug supplied with the equipment replaced by a technician with an approved plug suited for this power supply system.
-  The device is grounded when delivered. If the original plug is changed, you must make absolutely sure that the protective ground wire is connected to the new plug.

SAFETY INSTRUCTIONS

-  Please follow all safety and accident prevention regulations in effect in the laboratory.
-  Before plugging the device into the power supply circuit, ensure that the mains voltage matches the data on the rating plate.
-  Repairs may be made only by a Heidolph Instruments authorized technician.
-  Caution when used near easily inflammable and explosive materials. Please read the safety data sheets. The motors are designed to operate without sparking but the device is not protected from explosion.
-  Connect the disperser to a grounded mains socket only.
-  Make sure the device is placed on a stable surface.
-  When operating, wear suitable clothing and eye protection.
-  Caution: The tool and fluid for dispersion are heated by the energy input.


ASSEMBLY

Caution: Never operate the SilentCrusher M with dispersion tool when the dispersion tool is not immersed in fluid. Operating the device dry will destroy the PTFE bearing (AA, see Chapter 0) and the rotor (BB, see Chapter 0).

1. Set up Stand (if ordered as an option)

Screw the stand rod (A) into the stand plate (B). Using the wrench (C) (included) tighten stand rod (A) on the level areas (D).

Place the sleeve coupling (E) over the stand rod (A) and tighten it with set screw 1 (F).

2. Mount Drive Unit

Insert the retainer rod (G) for the drive unit into the hole in the upper sleeve coupling (E) and tighten with the set screw 2 (H).

3. Mount Vessel with Vessel Retainer (if ordered as an option)

The vessel retainer comes with 2 retainer rings (I and J). Select the right size for your vessel.


Size	For Vessel Diameter [mm]
I (small)	20 - 60
J (large)	60 - 105

Insert the vessel retainer (I or J) with the retainer rod (G) into the hole in the sleeve coupling (E) and tighten with the set screw 2 (H). Loosen the knurled screw 3 (L) and pull the retainer ring (I or J) forward. Insert the vessel (M) from above, or if the shape of the vessel requires it, from the bottom.

Then press the retainer ring downward as shown (support the sleeve coupling (E) with your index finger while doing this) and tighten the knurled screw (L). The vessel is now mounted.

Check that the vessel is mounted securely before operating.


4. Insert Dispersion Tool (Standard)

Insert the dispersion tool (N) into the hole (O) in the drive unit up to the stop and secure it with the set screw (P).

Note: At the upper end of the rotor shaft, the dispersion tools have a strong magnet that is responsible for inducing rotation.


Caution when doing manipulations in the vicinity of iron parts or other dispersion tools. The parts will attract each other forcefully.


5. Insert Sealed Dispersion Tool

 As delivered, the drive unit is equipped for the standard dispersion tool (non-sealed). To use the sealed dispersion tool, the drive unit must be converted.

To convert the drive unit, place it in a horizontal position. Then insert the screwdriver (Q) (included) into the slit in the inner ring of the ball bearing 1 (R) and unscrew it clockwise. Remove the ball bearing 1 (R) for use with the standard dispersion tool. Then return the drive unit to an upright position.

Insert the sealed dispersion tool (N) into the hole (O) in the drive unit (U) up to the stop and then screw it on by rotating it counterclockwise. Secure with set screw (P).

 The following rules apply to the conversion of the drive unit:

- For non-sealed dispersion tools, insert ball bearing 1 (R).
- For sealed dispersion tools, remove ball bearing 1 (R).

6. Positioning the Dispersion Tool in the Fluid

The dispersion tool, with the drive unit and the vessel retainer, can be positioned correctly by loosening the knurled screws 1 and 2 (F and H) on the sleeve coupling and moving the retainer rod. Retighten the knurled screws 1 and 2 (F and H) on the sleeve coupling (E).

7. Operating under Pressure or in a Vacuum with the Sealed Dispersion Tool


When operating under pressure or in a vacuum with the sealed tool, always position the stator shaft in the fluid so that the upper hole (TT) in the stator shaft is always above the level of the fluid.

If the hole (TT) is below the level of the fluid, fluid can enter the upper area of the dispersion tool due to the change in pressure conditions and damage it.

8. Connecting the Power Supply (unit) (S)

Connect the Power Supply (unit) (S) to the drive unit (U) with the small plug (T).


Connect the power connection cable (V) to the receptacle (X) with the plug (W). Place the switch (Y) in the “0” position and turn the knob (Z) all the way to the left stop.

9. Start up Device

Caution: Never operate the SilentCrusher M with the dispersion tool when the dispersion tool is not immersed in fluid. Operating the device dry will destroy the PTFE bearing (AA, see Chapter 0) and the rotor (BB, see Chapter 0).

Connect the power connection cable (V) to the power supply. The green control light (CC) on the power supply (unit) (S) will light. Turn the device on with the switch (Y). The display (DD) will light.

Be sure the vessel is filled with the sample for dispersion and that the dispersion tool is immersed at least 25 mm in the fluid.

Press the Start-Stop key (EE) to start the dispersion process. The LED (FF) will glow green.

Set the desired dispersion speed on rotating knob (Z). The display (DD) will show the speed (The display takes place in steps of 500).

If the dispersion tool is not in the drive unit, the speed indicator will remain at 0.0.

Caution: Due to its low power consumption, the device does not have a mains switch. If the device is connected to the mains voltage (standby), the device consumes only 6 Watts.

To disconnect the device from the power supply system, unplug the mains plug.

When operating under pressure or in a vacuum, always position the stator shaft with the sealed tool in the fluid so that the upper hole (TT) in the stator shaft is always above the level of the fluid.

If the hole (TT) is below the level of the fluid, it can enter the upper range of the dispersion tool and damage it.

Please note minimum and maximum immersion depth!

* minimal immersion depth using vacuum

	Min.	Min. (vac)*	Max.
8DF	25	35	50
12DF	25	35	75
12DG	25	35	75
18DF	25	35	135
18DG	25	35	135

	Min.	Max.
6F	20	50
8F	25	55
12F	25	80
12G	25	80
18F	25	140
18G	25	140
22F	30	140
22G	30	140

9.1 Speed limitation of the drive

The tools 22 F and 22 G can be operated with a maximum speed of 20,000 rpm. If this speed is exceeded, the drive automatically limits the speed after a brief period.

The limitation to 20,000 rpm is portrayed by the blinking right decimal point on the display (DD) and can be reset by reducing the desired speed <20,000 rpm.

With an appropriate load, speed limitation can also activate with tools <22 F/G.


10. Tool Replacement

Always disconnect the device using the switch (Y) before replacing the tool.

Standard Dispersion Tool (1)

Loosen the set screw (P) and pull the dispersion tool (N) downward. Insert a new dispersion tool and re-secure with the set screw (P). An elastomer O ring (QQ) inserted in opening (O) provides a secure fit, without play, for the dispersion tool.

Sealed Dispersion Tool (2)

Loosen the set screw (P) and screw the dispersion tool (N) counterclockwise out of the opening and pull downward. Insert a new tool up to the stop and screw in clockwise. Then re-secure with the set screw (P).

 **If you switch from a sealed to an unsealed dispersion tool or vice versa, the drive unit must be converted. See Installation, Chapter 5.**

11. Length of Operation

The drive unit and dispersion tool are not designed for constant operation. In particular, the PTFE bearing (AA) could be damaged by the high rotational speeds of 26000 rpm.

Keep dispersing time at a minimum, to save the PTFE bearing.

The maximum operation time is dependent on the heating-up of the drive (as consequence of load on the tool).

Operation time after first dispersion process is a minimum of 30 minutes. For subsequent dispersion processes, operation time can be shortened according to the heating-up of the device. Longer pauses in operation lead to longer operation times.

With sealed dispersion tools slightly reduced operation times have to be considered.

Times reported here are related to water as dispersing liquid at an ambient temperature of 22°C.

The electronics and coil in the drive unit are also monitored by temperature sensors. If the permissible temperature threshold is exceeded, the drive switches off. The LED (GG) in the service part of the drive unit will glow red. After cooling (about 15 min.), this LED will go out and dispersion can be continued by pressing the Start-Stop key (EE).

When the PTFE bearing in the dispersion tool becomes worn, it must be replaced. See the Cleaning and Maintenance chapter.

Wear can be recognized from the increased play between the rotor (BB) and the PTFE bearing (AA) or the increased level of noise.

The dispersion tools will heat in the area of the interlock casing (JJ) (see Chap. 0); in particular, the sealed dispersion tools can heat to the point that there is a risk of combustion. Take precautions when replacing the tool and dismantling the tool.

CLEANING AND MAINTENANCE

12. Drive Unit

For **cleaning**, you can wipe down the housing and surface of the device with a moist cloth (mild soap solution).

Note:

Never use chlorine bleach, chlorine-forming cleansers, abrasives, ammonia, steel wool or cleaning agents with metal particles. They can damage the surface of the device.

The device is maintenance-free. Any repairs required must be performed by a Heidolph Instruments authorized technician. Contact your Heidolph Instruments dealer or a representative of Heidolph Instruments (see page 54).

13. Dispersion Tool


To avoid contamination, the dispersion tool must be cleaned after every dispersing process.

13.2 Cleaning:


Clean the dispersion tool with a suitable cleaning fluid (depending on the medium processed). Run the dispersion tool in the cleaning fluid so that the internal parts are also cleaned. If the tool is extremely dirty, the dispersion tool can be disassembled for cleaning.

13.3 Disassembly, Standard Dispersion Tool:

Insert the auxiliary tool (HH) into the cutting area of the tool in the manner shown. Note: the auxiliary tool (HH) has three different size insertion areas. These are designed for the 3 sizes of tools: 8 / 12 / 18. Select the right size for your tool.

Then unscrew the cutting rotor (II) from the rotor shaft (BB) clockwise. Remove the cutting rotor (II) and the rotor shaft (BB). To remove the PTFE bearing (AA) insert the cutting rotor on other side of rotor shaft (BB) and gently push bearing out. The elastomer O ring (PP) on the rotor shaft provides a play-free seat for the rotor shaft on the ball bearing 1 (R).


Dispersion Tool 6F

Disassembly Dispersion Tool 6F

Use allen key (HH) to loosen screw (DD). Remove the rotor shaft (BB). To remove PTFE bearing (AA) stick rotor shaft into opposite side of stator shaft (MM) and push out bearing.

The elastomer o-ring (PP) on the rotor shaft provides a play-free seat for the rotor shaft on the ball bearing (R)

Assembly of dispersion tool 6F

Place the PTFE-bearing (AA) on the rotor shaft (BB). Then place the rotor shaft (BB) plus the PTFE bearing (AA) on the stator shaft (MM). Use screw (DD) to fix rotor and secure with allen key.


13.4 Disassembly, Sealed Dispersion Tool:


The dispersion tools will heat in the area of the interlock casing (JJ); in particular, the sealed dispersion tools can heat to the point that there is a risk of combustion. Take precautions when replacing the tool and dismantling the tool.

Unscrew the interlock casing (JJ) using the wrenches (RR) and (SS) (included) counterclockwise. Use the wrench surfaces (NN) provided on both parts for this.

Pay attention to the elastomer O ring (OO) inserted to provide a seal in the interlock casing (JJ);


it must not be damaged and is responsible for sealing the dispersion tool.

Insert the auxiliary tool (HH) into the cutting area of the dispersion tool in the manner shown.

Caution: The auxiliary tool (HH) has three different size insertion areas. These are designed for the 3 sizes of tools: 8 / 12 / 18. Select the right size for your tool. Then unscrew the cutting rotor (II) from the rotor shaft (BB) clockwise. Remove the cutting rotor (II) and the rotor shaft (BB). To remove the PTFE bearing (AA) insert the cutting rotor on other side of rotor shaft (BB) and gently push bearing out.


The ball bearing 2 (KK) can be removed from the interlock casing (JJ) for cleaning. Insert the screwdriver (Q) (included) into the slot in the inner ring of the ball bearing 2 (KK) and unscrew it clockwise.

13.5 Sterilization of the Dispersion Tool:

This dispersion tool can be sterilized as follows:
 By steam sterilization or wet chemicals at temperature up to 130°C; dry up to 180°C

Caution: Check the ball bearing 2 (KK) of the sealed dispersion tool after cleaning or sterilization: The balls of the ball bearing must always be lubricated and never be allowed to run dry since this could damage them immediately. Re-lubricate the ball bearings 2 (KK) after every cleaning or sterilization with the special lubricant (LL) (included).

13.6 Assembly of the Dispersion Tool:

Place the PTFE bearing (AA) on the rotor shaft (BB). Then place the rotor shaft (BB) plus the PTFE bearing (AA) on the stator shaft (MM). Place the cutting rotor (II) on the other end of the stator shaft (MM), with the threading first. Insert the auxiliary tool (HH) into the cutting rotor (II) and tighten it by rotating the rotor shaft (BB) clockwise.

For sealed dispersion tools only: Screw the ball bearing 2 (KK) back onto the interlock casing (JJ). Screw the interlock casing (JJ) onto the stator shaft (MM) and tighten with the wrench.

SPARE AND WEAR PARTS

Number	Description	Part Number
QQ	O Ring - Opening Drive Unit	23-08-06-03-17
R	Ball bearing 1 complete	11-001-595-13
Q	Screwdriver	04-06-01-21-13

Dispersion Tool 6 F

PP	o-ring rotor shaft (Viton)	23-08-06-03-12
BB	rotor 6F	22-02-06-01-84

AA	PTFE-bearing 6	25-05-01-03-07
MM	stator 6F	22-02-06-01-85
DD	Screw M3x4	02-01-19-03-15
HH	Allen key 2,5mm	02-07-01-01-04

Dispersion Tool 8 F

AA	PTFE Bearing	25-05-01-04-03
MM	Stator Shaft	22-02-06-01-50
BB	Rotor Shaft	22-02-06-01-60
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-63
HH	Auxiliary Tool	04-06-01-21-12

Dispersion Tool 8 DF

AA	PTFE Bearing	25-05-01-04-03
MM	Stator Shaft	22-02-06-01-51
BB	Rotor Shaft	22-02-06-01-60
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-63
JJ	Interlock Casing	22-02-06-01-66
OO	O Ring Interlock Casing (Viton)	23-08-06-03-16
KK	Ball Bearing, 2 complete	11-001-595-14
Q	Screwdriver	04-06-01-21-13
HH	Auxiliary Tool	04-06-01-21-12
LL	Special Lubricant	11-001-595-15
RR	Wrench 27	02-07-04-01-03
SS	Wrench 19	02-07-04-01-04

Dispersion Tool 12 F

AA	PTFE Bearing	25-05-01-05-02
MM	Stator Shaft	22-02-06-01-52
BB	Rotor Shaft	22-02-06-01-61
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-64
HH	Auxiliary Tool	04-06-01-21-12

Dispersion Tool 12 G

AA	PTFE Bearing	25-05-01-05-02
MM	Stator Shaft	22-02-06-01-54
BB	Rotor Shaft	22-02-06-01-61
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-64
HH	Auxiliary Tool	04-06-01-21-12

Dispersion Tool 12 DF

AA	PTFE Bearing	25-05-01-05-02
MM	Stator Shaft	22-02-06-01-53
BB	Rotor Shaft	22-02-06-01-61
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-64
JJ	Interlock Casing	22-02-06-01-66
OO	O Ring Interlock Casing (Viton)	23-08-06-03-16
KK	Ball Bearing, 2 complete	11-001-595-14
Q	Screwdriver	04-06-01-21-13
HH	Auxiliary Tool	04-06-01-21-12
LL	Special Lubricant	11-001-595-15
RR	Wrench 27	02-07-04-01-03
SS	Wrench 19	02-07-04-01-04

Dispersion Tool 12 DG

AA	PTFE Bearing	25-05-01-05-02
MM	Stator Shaft	22-02-06-01-55
BB	Rotor Shaft	22-02-06-01-61
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-64
JJ	Interlock Casing	22-02-06-01-66
OO	O Ring Interlock Casing (Viton)	23-08-06-03-16
KK	Ball Bearing, 2 complete	11-001-595-14
Q	Screwdriver	04-06-01-21-13
HH	Auxiliary Tool	04-06-01-21-12
LL	Special Lubricant	11-001-595-15
RR	Wrench 27	02-07-04-01-03
SS	Wrench 19	02-07-04-01-04

Dispersion Tool 18 F

AA	PTFE Bearing	25-05-01-06-01
MM	Stator Shaft	22-02-06-01-56
BB	Rotor Shaft	22-02-06-01-62
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-65
HH	Auxiliary Tool	04-06-01-21-12

Dispersion Tool 18 G

AA	PTFE Bearing	25-05-01-06-01
MM	Stator Shaft	22-02-06-01-58

BB	Rotor Shaft	22-02-06-01-62
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-65
HH	Auxiliary Tool	04-06-01-21-12

Dispersion Tool 18 DF

AA	PTFE Bearing	25-05-01-06-01
MM	Stator Shaft	22-02-06-01-57
BB	Rotor Shaft	22-02-06-01-62
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-73
JJ	Interlock Casing	22-02-06-01-66
OO	O Ring Interlock Casing (Viton)	23-08-06-03-16
KK	Ball Bearing, 2 complete	11-001-595-14
Q	Screwdriver	04-06-01-21-13
HH	Auxiliary Tool	04-06-01-21-12
LL	Special Lubricant	11-001-595-15
RR	Wrench 27	02-07-04-01-03
SS	Wrench 19	02-07-04-01-04

Dispersion Tool 18 DG

AA	PTFE Bearing	25-05-01-06-01
MM	Stator Shaft	22-02-06-01-59
BB	Rotor Shaft	22-02-06-01-62
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-73
JJ	Interlock Casing	22-02-06-01-66
OO	O Ring Interlock Casing (Viton)	23-08-06-03-16
KK	Ball Bearing, 2 complete	11-001-595-14
Q	Screwdriver	04-06-01-21-13
HH	Auxiliary Tool	04-06-01-21-12
LL	Special Lubricant	11-001-595-15
RR	Wrench 27	02-07-04-01-03
SS	Wrench 19	02-07-04-01-04

Dispersion Tool 22 F

AA	PTFE Bearing	25-05-01-06-01
MM	Stator Shaft	22-02-06-01-74
BB	Rotor Shaft	22-02-06-01-62
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-72
HH	Auxiliary Tool	04-06-01-21-12

Dispersion Tool 22 G

AA	PTFE Bearing	25-05-01-06-01
MM	Stator Shaft	22-02-06-01-75
BB	Rotor Shaft	22-02-06-01-62
PP	O Ring Rotor Shaft (Viton)	23-08-06-03-12
II	Cutting Rotor	22-02-06-01-72
HH	Auxiliary Tool	04-06-01-21-12

DISASSEMBLY, TRANSPORT AND STORAGE

Disassembly

Remove the mains plug and disconnect the connecting socket from the Power Supply to the drive unit.

Remove the dispersion tool from the drive unit.

Remove the drive unit from the stand; if necessary, remove the vessel retainer from the stand as well.

Transport and Storage

1. The device and its parts should be stored preferably in the original packaging or in another suitable container, to prevent damage during transport. The best material to seal the packing is adhesive tape.
2. Keep the device in a dry location.

Warning

Avoid shaking and shocks to the device during transport.

DISPOSAL

Dispose of used or defective devices properly in a collection location.

Separate recyclable material into metal, glass, plastic, etc.

Dispose of the old packing material in an environmentally friendly manner (material separation).

TROUBLESHOOTING

- Device does not disperse
 - Check electrical connections (power supply system and power supply)
 - Dispersion tool not inserted into drive unit up to stop
 - Maximum temperature reached in drive unit, red LED goes on (see Service and Operation)
 - Rotor not mounted in tool
 - Rotor sluggish (check PTFE bearing and ball bearing, see Cleaning and Maintenance)
- Loud noise when dispersing
 - PTFE bearing worn, must be replaced (see Cleaning and Maintenance, Chapters 13.2 and 13.5)
 - Ball bearing in drive unit worn, must be replaced (see Assembly, Chapter 5)
 - Ball bearing in sealed dispersion tool worn, must be replaced (see Cleaning and Maintenance, Chapters 13.3 and 13.5)
- Attrition (PTFE particles) at bearing
 - Tools were operated without or with too little liquid; clean tool, if necessary exchange bearing. Please note the immersion depths (p. 39).

If a malfunction occurs that cannot be corrected with the above instructions, please inform your authorized Heidolph Instruments dealer immediately.

TECHNICAL DATA

SilentCrusher M

Characteristics	Values
Connection Voltage	100 – 230 V 50 / 60 Hz
Power Consumption (Rated Output)	130 W
RPMs	5000 – 26000 rpm, referred to water
Operating Modes	Short-term operation
Drive Protection	Overheating protection switches drive unit off in case of overheating; red control light indicates this condition. Reset with start-stop key.
Control Light for “Device On”	CONTROL LIGHT (green) on power supply
Dimensions (in mm) Drive Unit	Width 68; Height 160; Depth 75
Retainer Rod	10 x 100
Dimensions (in mm) Power Supply	Width 145; Height 43; Depth 75
Total weight	1.9 kg
Operating Temperature	0 – 40 °C
Storage Temperature	- 20 to + 80 °C
Relative Humidity	85 % (without condensation)
Type of Protection according to EN 60529	IP 30
Level of Protection according to VDE 0100	I

Dispersion Tool 6 F

Characteristics	Values
Part Number	596-06010-00
Stator Diameter [mm]	5,3
Rotor Diameter [mm]	3,2
Gap Width Generator/Rotor [mm]	0,3
Width of Slot [mm]	0,6
Total Shaft Length [mm]	60
Maximum Immersion Depth [mm]	48
Minimum Immersion Depth [mm]	20
Volume [ml]	0,8-30
Peripheral Speed [m/sec]	4,4
Initial Particle Size [mm]	< 0,4
Material	PTFE / 1.4435
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 8 F

Characteristics	Values
Part Number	596-08010-00
Stator Diameter [mm]	8.2
Rotor Diameter [mm]	6.6
Gap Width Generator/Rotor [mm]	0.3
Width of Slot [mm]	0.6
Total Shaft Length [mm]	80
Maximum Immersion Depth [mm]	55
Minimum Immersion Depth [mm]	25
Volume [ml]	3 -50
Peripheral Speed [m/sec]	9
Initial Particle Size [mm]	< 0.4
Material	PTFE / 1.4435
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 8 DF

Characteristics	Values
Part Number	596-08110-00
Stator Diameter [mm]	8.2
Rotor Diameter [mm]	6.6
Gap Width Generator/Rotor [mm]	0.3
Width of Slot [mm]	0.6
Total Shaft Length [mm]	75
Maximum Immersion Depth [mm]	50
Minimum Immersion Depth [mm]	25
Minimum Immersion Depth [mm] under vacuum	35
Volume [ml]	3 - 50
Peripheral Speed [m/sec]	9
Initial Particle Size [mm]	< 0.4
Material	PTFE / 1.4435
Maximum Permissible Pressure	6 bar
Permissible Vacuum	20 mbar
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 12 F

Characteristics	Values
Part Number	596-12010-00
Stator Diameter [mm]	12.0
Rotor Diameter [mm]	9.0
Gap Width Generator/Rotor [mm]	0.3
Width of Slot [mm]	0.6
Total Shaft Length [mm]	120
Maximum Immersion Depth [mm]	80
Minimum Immersion Depth [mm]	25
Volume [ml]	25 - 250

Peripheral Speed [m/sec]	12.2
Initial Particle Size [mm]	< 2.0
Material	PTFE / 1.4435
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 12 G

Characteristics	Values
Part Number	596-12020-00
Stator Diameter [mm]	12.0
Rotor Diameter [mm]	9.0
Gap Width Generator/Rotor [mm]	0.3
Width of Slot [mm]	2.0
Total Shaft Length [mm]	120
Maximum Immersion Depth [mm]	80
Minimum Immersion Depth [mm]	25
Volume [ml]	3 - 50
Peripheral Speed [m/sec]	8
Initial Particle Size [mm]	< 0.4
Material	PTFE / 1.4435
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 12 DF

Characteristics	Values
Part Number	596-12110-00
Stator Diameter [mm]	12.0
Rotor Diameter [mm]	9.0
Gap Width Generator/Rotor [mm]	0.3
Width of Slot [mm]	0.6
Total Shaft Length [mm]	115
Maximum Immersion Depth [mm]	75
Minimum Immersion Depth [mm]	25
Minimum Immersion Depth [mm] under vacuum	35
Volume [ml]	25 - 250
Peripheral Speed [m/sec]	12.2
Initial Particle Size [mm]	< 2.0
Material	PTFE / 1.4435
Maximum Permissible Pressure	6 bar
Permissible Vacuum	20 mbar
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 12 DG

Characteristics	Values
Part Number	596-12120-00
Stator Diameter [mm]	12.0
Rotor Diameter [mm]	9.0
Gap Width Generator/Rotor [mm]	0.3
Width of Slot [mm]	2.0
Total Shaft Length [mm]	115
Maximum Immersion Depth [mm]	75
Minimum Immersion Depth [mm]	25
Minimum Immersion Depth [mm] under vacuum	35
Volume [ml]	25 - 250
Peripheral Speed [m/sec]	12.2
Initial Particle Size [mm]	< 2.0
Material	PTFE / 1.4435
Maximum Permissible Pressure	6 bar
Permissible Vacuum	20 mbar
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 18 F

Characteristics	Values
Part Number	596-18010-00
Stator Diameter [mm]	18.0
Rotor Diameter [mm]	12.5
Gap Width Generator/Rotor [mm]	0.4
Width of Slot [mm]	1.5
Total Shaft Length [mm]	180
Maximum Immersion Depth [mm]	140
Minimum Immersion Depth [mm]	25
Volume [ml]	100 - 1500
Peripheral Speed [m/sec]	17
Initial Particle Size [mm]	< 4.0
Material	PTFE / 1.4435
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 18 G

Characteristics	Values
Part Number	596-18020-00
Stator Diameter [mm]	18.0
Rotor Diameter [mm]	12.5
Gap Width Generator/Rotor [mm]	0.4
Width of Slot [mm]	3.0
Total Shaft Length [mm]	180
Maximum Immersion Depth [mm]	140
Minimum Immersion Depth [mm]	25
Volume [ml]	100 - 1500
Peripheral Speed [m/sec]	17
Initial Particle Size [mm]	< 4.0
Material	PTFE / 1.4435
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 18 DF

Characteristics	Values
Part Number	596-18110-00
Stator Diameter [mm]	18.0
Rotor Diameter [mm]	12.5
Gap Width Generator/Rotor [mm]	0.4
Width of Slot [mm]	1.5
Total Shaft Length [mm]	175
Maximum Immersion Depth [mm]	135
Minimum Immersion Depth [mm]	25
Minimum Immersion Depth [mm] under vacuum	35
Volume [ml]	100 - 1500
Peripheral Speed [m/sec]	17
Initial Particle Size [mm]	< 4.0
Material	PTFE / 1.4435
Maximum Permissible Pressure	6 bar
Permissible Vacuum	20 mbar
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 18 DG

Characteristics	Values
Part Number	596-18120-00
Stator Diameter [mm]	18.0
Rotor Diameter [mm]	12.5
Gap Width Generator/Rotor [mm]	0.4
Width of Slot [mm]	3.0
Total Shaft Length [mm]	175
Maximum Immersion Depth [mm]	135
Minimum Immersion Depth [mm]	25
Minimum Immersion Depth [mm] under vacuum	35
Volume [ml]	100 - 1500
Peripheral Speed [m/sec]	17
Initial Particle Size [mm]	< 4.0
Material	PTFE / 1.4435
Maximum Permissible Pressure	6 bar
Permissible Vacuum	20 mbar
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 22 F (max. speed 20.000 rpm)

Characteristics	Values
Part Number	596-22010-00
Stator Diameter [mm]	22.0
Rotor Diameter [mm]	16.5
Gap Width Generator/Rotor [mm]	0.4
Width of Slot [mm]	1.5
Total Shaft Length [mm]	180
Maximum Immersion Depth [mm]	140
Minimum Immersion Depth [mm]	30
Volume [ml]	100 - 2000
Peripheral Speed [m/sec]	17.3 at 20.000 rpm
Initial Particle Size [mm]	< 4.0
Material	PTFE / 1.4435
Sterilization	Up to 130°C with wet chemicals; up to 180°C dry

Dispersion Tool 22 G (max. speed 20.000 rpm)

Characteristics	Values
Part Number	596-22020-00
Stator Diameter [mm]	22.0
Rotor Diameter [mm]	16.5
Gap Width Generator/Rotor [mm]	0.4
Width of Slot [mm]	3.0
Total Shaft Length [mm]	180
Maximum Immersion Depth [mm]	140
Minimum Immersion Depth [mm]	30
Volume [ml]	100 - 2000
Peripheral Speed [m/sec]	17.3 at 20,000 rpm
Initial Particle Size [mm]	< 4.0
Material	PTFE / 1.4435

Sterilization	Up to 130°C with wet chemicals; up to 180°C dry
----------------------	---

WARRANTY, LIABILITY AND COPYRIGHTS

Warranty

Heidolph Instruments warrants the products described here (except for wear parts) for three years from delivery from the manufacturer's warehouse. This warranty covers defects in materials and manufacturing.

Transportation damages are excluded.

For warranty claims, contact Heidolph Instruments (Tel.: (+49) 9122 - 9920-68) or your Heidolph Instruments dealer. If it involves a defect in materials or manufacturing, the device will be repaired or replaced under the warranty, free of charge.

Heidolph Instruments does not cover damage caused by improper handling.

Any change in this warranty must be confirmed by Heidolph Instruments in writing in each individual case.

Disclaimer of Liability

Heidolph Instruments does not accept any liability for improper handling and use of this device. Consequential damages are excluded.

Copyright

Heidolph Instruments holds the copyright to all images and text in this operating manual.

QUESTIONS / REPAIRS

If, after reading this operating manual, you still have **questions** about the installation, operation or maintenance of your device, please submit them to the following address.

For **repairs**, first call Heidolph Instruments directly (Tel.: (+49) 9122 - 9920-69) or your authorized Heidolph Instruments dealer.

Note:

Only after you have received approval, please send your device to:

Heidolph Instruments GmbH & Co. KG
Vertrieb Labortechnik
Walpersdorfer Str. 12
D-91126 Schwabach / Germany
Tel.: +49 – 9122 - 9920-74
Fax: +49 – 9122 - 9920-65
E-Mail: sales@heidolph.de

Safety Warning

When returning equipment for repair that has been in contact with hazardous materials, please indicate:

- The most accurate possible *material data* on the corresponding medium
- *Protective measures* required for safe handling by our receiving and maintenance personnel
- *Coding* on the packing pursuant to the Hazardous Materials Ordinance.

EC DECLARATION OF CONFORMITY

We declare that this product meets the following standards and normative documents:

Electromagnetic Compatibility Guidelines (89/336/EEC):

EN 61326: 1997 + A1:1998 + A2:2001
 EN 61000-3-2:2000
 EN 61000-3-3:1995
 EN 61326: 1997 + A1:1998 + A2:2001
 EN 61000-4-2:1995
 EN 61000-4-3:1996
 EN 61000-4-4:1995
 EN 61000-4-5:1995
 EN 61000-4-6:1996

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© Heidolph Instruments GmbH & Co. KG