Safe

Heating and Mixing

Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

Operating Manual

Page 2

Magnetic Stirrer MR Hei-End
Operating manual must be read before initial start-up.
Please follow the safety instructions provided.
Please keep for future reference.

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About this Document	4
Version	4
About this manual	4
Terms and signs used.....	5
Safety Instructions	6
General safety instructions.....	6
EU Declaration of Conformity.....	6
Intended use.....	6
Installation / Electrical safety	6
Personnel qualification.....	7
Operating company's obligations.....	7
Installation site	7
Changes to the unit	7
Personnel safety	8
Safety during use	8
Disposal	8
Device Description	9
Device overview	9
MR Hei-End overall view	9
MR Hei-End control panel.....	9
Start-Up	10
Setting up device	10
Connect/disconnect power cord	10
Switch unit on and off	11
Operation	12
Place vessels.....	12
Start mixing and heating	12
Set speed and start rotation.....	12
Set temperature	13
Start heating.....	15

Cleaning and Maintenance	16
Cleaning.....	16
Maintenance	16
Troubleshooting	17
Possible errors and how to resolve	17
Assembly	19
Electrical connections	19
Mounting periphery devices	20
Connect temperature sensor Pt 1000.....	20
Connect PC	21
Configure Interfaces.....	22
Use interface analog 0-10V	23
Disassembly and Storage	24
Dismantling, transportation and storage	24
Dismantling	24
Transportation and storage	24
Accessories and Spare Parts	25
Scope of delivery.....	25
Accessories.....	25
Attachments	26
Technical data	26
Service	27
Contact / Technical Service	27
Warranty	28
Confirmation of condition of unit.....	29

➤ Version

Version

This manual describes the function, operation and maintenance of the magnetic stirrer MR Hei-End

Version	Alteration date
1.0	11.2017

➤ About this manual

This operating manual is a component of the device described here.

- Please read this manual carefully and obey all safety and warning notices.
- Ensure that every operator reads this manual.
- Ensure that this manual is accessible for every operator.
- Pass on the operating manual to the subsequent owner.

For the current version of this operating manual in pdf format, please go to:

<http://www.heidolph-instruments.com/Support/Operation Manuals/Magnetic Stirrers/All MR Hei- models/MR Hei-End>

In addition please observe the regional regulations.

➤ Terms and signs used

In this manual you will find standardized terms and signs to warn you of possible dangers or give important hints. Please take special notice to these signs to avoid personal injury or damage to the unit.

The following terms and type of signs (basic drawing elements) are used:

Symbol	Additional term / Description
Warning signs 	The yellow triangle indicates hazardous situations. It is used in combination with the following terms: DANGER: Indicates a hazardous situation which, if not avoided, will result in irreversible injury or death. WARNING: Indicates a hazardous situation which, if not avoided, could result in serious injury or material damage. CAUTION: Indicates a hazardous situation which, if not avoided, may result in a minor injury or material damage.
Prohibitory signs 	Not to do: The red circle indicates a situation that should be avoided under all circumstances to avoid injury or damage.
Mandatory signs 	To do: The blue circle indicates important information. Please obey to avoid property damage.

Other signs used

Symbol	Description
→	Handling instruction, action required
✓	Result of action
▪ ▪ ▪	List of information
a. b. c.	List of variants

➤ General safety instructions

EU Declaration of Conformity

This device complies with the following EC-Directives*:

- 2006/42/EC Machinery Directive
- 2014/30/EU Electromagnetic Compatibility Directive

*Also see attachment „EU Declaration of Conformity“.

The device has been constructed according to state-of-the-art technology and recognized safety regulations. However, risks may still arise during installation, operation and maintenance.

→ Please ensure the operating manual is available at all times.

The device may only be used under the following circumstances:

- Only operate the device, if it is in full working order.
- Ensure all operators of the device possess the necessary safety and risk awareness.
- Operate the device according to the instructions stipulated in this manual only.
- If there is something you do not understand, or certain information is missing, ask your manager or contact the manufacturer.
- Do not do anything on the device without authorization.
- Only use the device according to its intended use.

Intended use

The device is intended for use by trained and authorized personnel only.

The device is suitable for the following use:	<ul style="list-style-type: none"> • heating • mixing
---	---

The operation for the device is suitable in one of the following locations:	in any research application
---	-----------------------------

Improper use

Any use which deviates from the device's intended use is considered to be improper. The manufacturer does not accept liability for any damages resulting from non-permitted uses. The risk is carried by the operator alone.

Installation / Electrical safety

- The device may only be connected when the mains voltage corresponds to the information on the type plate of the unit.
- The mains connection must be accessible at all times.
- Repairs may only be performed by a qualified electrician.
- Never operate the unit with a damaged power cord.
- Always turn the unit OFF and disconnect mains power before performing any maintenance or service.

Personnel qualification

- The device may only be operated by qualified persons.
- The device may only be operated by individuals who have been instructed in its proper use by qualified persons.
- The device may only be operated and maintained by persons who are of legal age.
- Other personnel may only work with the unit under continuous supervision of an experienced qualified operator.
- This manual must be read and understood by all persons working with the device.
- The personnel must have received special safety instructions in order to guarantee responsible and safe work procedures.

➤ Operating company's obligations

Installation site

- The device must be positioned in a suitable location.
- The device must be installed on a firm and stable surface.
- All screw connections must be securely tightened.
- Do not operate the device in proximity of highly inflammable and explosive substances.
- Installation and operation of the device is only permitted in facilities which are equipped with appropriate laboratory equipment (e.g. with extractor hoods).
- Only operate the device in conjunction with an extractor hood (at least 10-fold air change, with error monitoring), if working with potential harmful media (see DIN EN 14175 and DIN 12924).
- The device may only be operated in enclosed spaces and under the following environmental conditions:

Ambient temperature	5 - 31 °C at 80 % relative humidity 32 - 40 °C decrease linearly till max. 50 % relative humidity
Installation altitude	0 - 2,000 m (6,500 feet) height above sea level
Degree of pollution	2
Overvoltage category	II

- If the equipment is used in corrosive atmospheres, the service life of the equipment will decrease based on concentration, volumes, and frequency of exposures to these corrosive materials, for example concentrated Hydrochloric Acid (HCl).

Changes to the unit

- No unauthorized changes may be made to the unit.
- No parts may be used which have not been approved by the manufacturer.
- Unauthorized changes result in the EC Declaration of Conformity losing its validity, and the unit may no longer be operated.
- The manufacturer is not liable for any damage, danger or injuries that result from unauthorized changes or from operating the unit other than described in this manual.

Personnel safety

- Ensure that only qualified personnel operates the device.
- Observe the following regulations:
 - Laboratory guidelines
 - Accident prevention regulation
 - Ordinance on Hazardous Substances
 - Other generally accepted rules of safety engineering and occupational health
 - Local regulations

➤ Safety during use

- Beware of the effect of magnetic fields on cardiac pacemakers or magnetic sensitive data media. Observe the according warning in chapter "Start-up", "Setting up device".
- Wear the appropriate protective clothing when working on the device (protective glasses and, if necessary, safety gloves).
- Do not use the device in potentially explosive areas. The device is not protected against explosion. There is no explosion or ATEX protection available.
- Do not carry out work with naked flames in the vicinity of the device (risk of explosion).
- Do not operate or assemble devices in the vicinity which are emission or radiation sources (electromagnetic waves) for the frequency range ($3 \cdot 10^{11}$ Hz to $3 \cdot 10^{15}$ Hz).
- Do not operate or assemble appliances in the vicinity of the device which constitute emission or radiation sources for ionizing radiation or in the ultrasonic range.
- Do not operate the device where adiabatic compression or shock waves might occur (shock wave combustion).
- Do not use substances, where the energy input due to mixing might pose a hazard.
- Do not spill liquids over the device or any parts of it.
- Remove any accidentally spilled liquids immediately.
- Eliminate errors immediately.
- Always switch the device OFF after use.
- Avoid putting pressure on the display when you are not operating the device.
- Do not use abrasive material to clean the surface. Only wipe with damp cloth.

➤ Disposal

- Check the device components for hazardous substances and solvents.
- Clean all components before disposal.
- Dispose of the device according to the valid national legal regulations.
- Dispose of the packaging material in accordance with the valid national legal regulations.


➤ Device overview

MR Hei-End overall view

MR Hei-End control panel

➤ Setting up device

Warning: **Magnetic field!**

The magnetic field beneath the hotplate could cause irritations on your cardiac pacemaker, if you get too close.

Personnel with a cardiac pacemaker, implanted defibrillator or dosing pump should keep a minimum distance of 10 cm (4 inches) to the device.

Warning: **Risk of slipping device!**

As a result of vibrations generated during operation the device might slide off the table top.

The surface must be smooth, clean and temperature resistant and should be properly leveled.

- Locate the device on a stable, horizontal and temperature resistant surface.
- Clean surface and feet with a damp cloth and ethanol regularly.

➤ Connect/disconnect power cord

Connect power cord

The power cord comes with a three-wire plug with recess on the bottom side.

Unit plug

Power cord

The socket is located on the back side of the unit.

- ✓ The device is switched OFF.
- Connect unit plug to the unit socket with the recess facing down.
- Push plug up against socket until securely attached.
- Connect power plug to power socket.

Unit socket

Disconnect power cord

- Disconnect power plug from power socket first.
- Then disconnect plug from socket on unit.

➤ Switch unit on and off

The ON/OFF switch is situated at the right side of the device beneath the control panel.

Switch unit on

- Press the ON/OFF switch at the right side.
- ✓ The unit is switched on.

Switch unit off

- Press the ON/OFF switch at the left side.
- ✓ The unit powers off.

ON/OFF switch

➤ Place vessels

Caution:
Risk of breakage and spillage

If the vessels are not placed securely they may slip off and spill the sample.

Prior to start heating and mixing ensure that all vessels are securely placed.

Heating bath accessories for round flasks and beakers as well as attachments for water, gas and evaporating distributors are available.

- To place a number of vessels securely on plate use optional adaptors (see general catalogue).
- Single vessels with flat bottoms like beakers may be placed directly on hotplate.

A single flask should be arranged in the middle of the hotplate.
Several flasks should be distributed equally on the plate.

➤ Start mixing and heating

When using device inside conditioned cabinets, make reference to ambient conditions as stipulated in chapter "Technical data".

Set speed and start rotation

Warning:
Risk of poisoning!

Open vessels and too high speed may result in samples splashing.

Use vessels with narrow neck especially if sample is dangerous or toxic.

Adjust speed step by step until you have reached the required rpm settings.

Wear safety glasses and especially with dangerous and toxic samples appropriate protective clothing.


Operation

Speed can be set between 30 - 1,400 rpm.

- Place vessel with sample and stirring bar on hotplate.
- Turn item on with ON/OFF switch.
- Set speed with analog knob.
- ✓ Actual speed appears in the 3-digit digital display.
- ✓ Rotation starts immediately.

Temperature display: actual 36 °C

Speed display: actual 400 rpm

Set temperature

Warning:
Risk of burning!

The hotplate can reach temperatures up to 300 °C.

If you work with temperatures above 50 °C wear safety gloves and appropriate protective clothing.

Warning:
Risk of ignition!

Inflammable samples may ignite with high temperatures.

If you work with inflammable samples set temperature at least 25 °C lower than ignition point of sample.

There are two options to control heating:

- Without temperature sensor: temperature will be measured and controlled via hotplate.
- With temperature sensor: temperature will be measured and controlled via sensor directly in medium (see chapter "Assembly, Connect temperature sensor Pt 1000").

You can set following temperature range:

- Hotplate / temperature sensor: between 20 °C and 300 °C (lower degrees depending on ambient temperature).
- Safety delta: between 10 °C and 25 °C
If temperature raises above set temperature and exceeds safety delta, heating is automatically cut off.

Working without temperature sensor:


Due to differences in heat conductivity of various materials of vessels and medium the temperature of hotplate should be set 1½ to 4 times higher, than the desired medium temperature.

Operation

- ✓ Actual temperature is shown on top of display.
- Press button "Select" once.
- ✓ An arrow symbol appears at beginning of first line.
- ✓ Display changes from "hotplate temp" to "set hotplate".
- Set temperature in between 3 seconds using arrow buttons ▲ and ▼.
- Press button "Select" again in between 3 seconds.
- ✓ Arrow symbol appears at beginning of second line.
- Set delta for safety temperature in between 3 seconds using arrow buttons ▲ and ▼.

Set temperature: 100 °C

Delta safety temperature: Δt 25 °C

If no entries are made in between 3 seconds the last entry will be saved and the display switches to show actual temperatures.

Working with temperature sensor Pt 1000:


The display shows "set probe" instead of "set hotplate". Setting temperature is exactly the same procedure.


Operation

Start heating

- Press the heating button.
- ✓ The ring around heating symbol will illuminate green.
- ✓ The unit starts heating
- Set temperature in between 3 seconds using arrow buttons ▲ and ▼.


- ✓ A heating symbol ☰ will flash in display till set temperature is reached
- Press button "Select" to adjust heating temperature.
- Press heating button again to switch off heating.

The device is equipped with a residual heat warning.

- ✓ If device is switched off while temperature is still above 50 °C, the display will remain illuminated.
- ✓ A message will be displayed and actual temperature is shown.
- ✓ Display will switch off automatically when temperature lowers to 50° C and beneath.

A power cut will stop heating and mixing and put device on stand by. Set temperature and speed will be saved.

› Cleaning

All surfaces of the unit allow for cleaning with a damp cloth and if required a mild soap lotion.

Warning:
Risk of electric shock!

If the device is still connected to power during cleaning, you might get in contact with live parts due to accidentally entering liquids.

Prior to cleaning:

- ensure that the device is switched OFF
- and the plug is pulled from the electric socket.

Caution:
Damage to the device

Accidentally entering liquids will cause the device to fail .

Only wipe with a damp cloth. Never use a soaked cloth.

Attention:
Damage to the surfaces

All surfaces can be damaged by improper cleaning.

In any case do not use:

- Chlorine bleach or any cleaning agent including chlorine
- Ammonia
- Abrasive cleaning agents such as cleaning rags, scrubbing agents or any other agents which include metal components

As a result of vibrations generated during operation the device might slide off the table top.

To ensure a stable position clean surface and feet with a damp cloth and ethanol regularly.

› Maintenance

The unit is maintenance-free. Any necessary repair must be performed by an authorized Heidolph distributor.

Please contact Heidolph Instruments or your local Heidolph Instruments distributor.

➤ Possible errors and how to resolve

Error message	Possible reason	Troubleshooting
hot plate sensor break	Hotplate sensor broken	Contact your local Heidolph Instruments distributor
hot plate sensor shorted	Short circuit hot plate sensor	Contact your local Heidolph Instruments distributor
hot plate safety sensor break	Safety sensor broken	Contact your local Heidolph Instruments distributor
hot plate safety sensor shorted	Short circuit safety sensor	Contact your local Heidolph Instruments distributor
hot plate sensor difference >15°C	Hotplate sensor difference to high	Contact your local Heidolph Instruments distributor
hot plate sensor set safety ΔT	Set safety delta exceeded, indication of exothermal reaction	Increase safety temperature hot plate sensor "set safety ΔT "
external sensor break	External sensor broken	Replace external sensor
external sensor shorted	Short circuit external sensor	Replace external sensor
external safety sensor break	External safety sensor Sensor broken	Replace external sensor
external safety sensor shorted	Short circuit external safety sensor	Replace external sensor
external sensor difference >15°C	External sensor difference to high	Replace external sensor
external sensor set safety ΔT	Set safety delta exceeded, indication of exothermal reaction	Increase safety temperature external sensor "set safety ΔT "
triac IDxx	Triac error	Contact your local Heidolph Instruments distributor
motor	Motor defective, no rotation	Contact your local Heidolph Instruments distributor
communication fault IDxx	Communication fault between Triac and/or monitoring	Contact your local Heidolph Instruments distributor
Kernel error IDxx	Kernel error (operating system)	Contact your local Heidolph Instruments distributor
external probe outside media	External sensor not in contact with medium	Insert sensor at least 20 mm into medium

Troubleshooting

Malfunctions	Possible reason	Troubleshooting
Power LED does not enlighten	No power	<ul style="list-style-type: none"> ▪ Check size and shape of plug and compatibility to your electrical socket ▪ Check circuit breakers
	Light diode defective	Contact your local Heidolph Instruments distributor
No rotation function	No magnetic bars in vessel	Insert magnetic bars
	rpm set to "0"	Set speed
No heating up function	Set temperature exceeded	Let device cool down
Probe temperature lower than desired	Temperature controlled via hotplate	Raise temperature of hotplate

If you experience a malfunction which can not be resolved, please contact your authorized Heidolph Instruments distributor immediately.

Electrical connections

- Repairs of any kind are allowed by **qualified professional electricians only**. Any improper repair can result in a dangerous situation. Contact your local Heidolph Instruments distributor for any repair you may have.

Warning:
Electric shock!

This unit must only be connected to a grounded electrical socket.

The unit and voltage must match. The type label on the back side of the unit provides all voltage specifics.

When shipped, the unit is grounded. In case the original plug is changed, the new plug must have a protective conductor!

Rating plate

Color code for electrical connections:

Europe	
GREEN/ YELLOW	PE: Protective conductor (Earth)
BLUE	N: Neutral conductor
BROWN	P: Phase

North America	
GREEN	PE: Protective conductor (Earth)
WHITE	N: Neutral conductor
BLACK	P: Phase

The unit is connected with the power cord supplied to the electrical socket. It plugs into the housing on the rear side of the unit.

For countries utilizing a plug other than the standard one supplied:

- The plug supplied may only be changed by a professional electrician.
- If utilizing an adaptor ensure that it is approved by the local regulations.

➤ Mounting periphery devices

Connect temperature sensor Pt 1000

Caution:
Risk of cable damage

The cable of temperature sensor might get damaged by to high temperatures.

Ensure the cable does not get in contract with very hot materials. Please use clamping system with cable guide tube (optional accessori, see chapter "Accessories and Spare Parts").

- Switch off device.
- Turn locking ring on plug anti clockwise to the left.
- Insert plug into 4-pole socket on backside of device.
- Push locking ring towards device and turn clockwise, till it audible engages.

4-pole socket

Plug temperature sensor

For reliable measurement the temperature sensor has to be inserted at least 20 mm.

Connect PC

Use the RS232 interface and a suitable program (e.g. Microsoft® Hyperterminal), or a suitably programmed interface to enter and read out the parameters like temperature and speed.

Required connectors at interface cable:

- Device interface: SUB-HDD 15 connector.
- PC interface: SUB-D 9 connector.

(Also see chapter "Accessories and spare parts".)

Caution: Possible malfunction

A faulty connection can lead to malfunction of the interface.

Only connect those lines in accordance with drawing opposite.

Connection MR Hei-End

Connector SUB-HDD 15

e.g. McMurdo HDE15PTD

1 — TXD

2 — RXD

5 — GND

Connection PC (COM x)

Socket SUB-D 9

e.g. McMurdo SDE9ST

2 — RXD

3 — TXD

5 — GND

Connection at device

- Switch off the device.
- Remove protective cap and keep for future use.
- Connect SUB-HDD 15 to socket at rear of device.
- Connect SUB-D 9 to com interface at PC.
- ✓ When interface is active, "*REMOTE*" starts flashing in display.

Socket: Interface RS232

Configure Interfaces

- Configure interface to receive data in 7bit ASCII format.
- Each command string has to be finalized with CR (carriage return) and LF (line feed).

Interface parameters

9600 baud; parity: 1 even, 7 bits, 1 stop bit

Read command	Description	Output
IN_PV_1	Query actual temperature: probe temperature in °C	
IN_PV_2	Query actual temperature: probe safety sensor in °C	
IN_PV_3	Query actual temperature: hotplate in °C	
IN_PV_4	Query actual temperature: hotplate safety sensor in °C	
IN_PV_5	Query actual motor speed in rpm	
IN_SP_1	Query set temperature probe / hotplate in °C	
IN_SP_2	Query set temperature safety delta in °C	
IN_SP_3	Query set motor speed in rpm	
IN_MODE_2	Behavior after power cut	

Write Command	Description	Output
OUT_SP_1 y	Enter set temperature probe / hotplate in °C	
OUT_SP_2 y	Enter set temperature safety delta in °C	
OUT_SP_3 y	Enter set motor speed in rpm	
OUT_MODE_2 y	<ul style="list-style-type: none"> ▪ 0 = all OFF after power return ▪ 1 = heating / stirring ON after power return (set values remain valid) 	<ul style="list-style-type: none"> ▪ "0" ▪ "*"

Command	General commands: Description	Output
START_1(*)	Start heater control, activate remote – alternating display *Remote* (see footnote)	
START_2(*)	Start motor, activate remote – alternating display *Remote* (see footnote)	
STOP_1	Stop heater control	
STOP_2	Motor speed=0 Motor-stop	
RESET	Reset to initial state: heater off, motor off, exit remote mode	
STATUS	Request status	<ul style="list-style-type: none"> ▪ 0 = ManualMode ▪ 1 = RemoteMode Start (START_1 OR START_2) ▪ 2 = RemoteMode Stop (STOP_1 AND STOP_2) ▪ <0 Error code in case of fault

* Attention: control panel of device is disabled in remote mode

To activate commands

- command string for speed always has to be followed by "START_2"
- command string for heating has to be followed by "START_1" only the first time.

As soon as heating button is activated all further adjustments of temperature will be executed immediately.

Examples: remote operation

Description	Command string	Result
Heater to be adjusted to 150 °C	OUT_SP_1 150	Heating adjustment to 150 °C, only if: START_1 is entered
Next step: Heater to be adjusted to 250 °C	OUT_SP_1 250	Heating adjustment to 250 °C immediately after input
Motor to be adjusted to 1000 rpm	OUT_SP_3 1000	Motor adjustment to 1000 rpm, only if: START_2 is entered
Next step: Motor to be adjusted to 500 rpm	OUT_SP_3 500	Motor adjustment to 500 rpm, only if: START_2 is entered

Use interface analog 0-10V

Interpreted analog the interface can be used to output the actual medium temperature.

Pin assignment SUB-HDD plug as shown opposite

MR Hei-End

Plug SUB-HDD 15

e.g. McMurdo HDE15PTD

9 GND

7 Medium temperature output 0 - 10V

Output values:

- 0 °C = 0 V
- per 100 °C = 3,5 V

➤ Dismantling, transportation and storage

Dismantling

Warning: Risk of injury!

While the unit is plugged in you might be at harm if the device is accidentally switched to ON mode.

Prior to dismantling unit ensure that the unit is switched off and the plug is pulled from the electric socket.

- ✓ The unit is switched off.
- Disconnect plug from power socket.
- Disconnect plug from socket on the unit.
- Remove all vessels.
- Disassemble attachments.

Transportation and storage

- Store and transport the unit and its components only if they were emptied and cleaned,
- Store and transport the unit and parts of unit in the original packing material or alternatively in an adequate container to prevent damages.
- Seal the packing carefully against unauthorized or mistaken opening.
- Store the unit in a dry and frost-free place.

Improper storage and transportation may result in damages to the system and the mechanical components.

Avoid any kind of shocks during transportation.

➤ Scope of delivery

Component	Quantity	Product number
MR Hei-End (230 V)	1	505-50000-00
Operating Manual English / German	1	01-005-004-59
Warranty registration / Confirmation of condition	1	01-006-002-78
EC Declaration of Conformity	1	01-001-025-02
Power cord	1	depending on country

➤ Accessories

Component	Quantity	Product number
Temperature sensor Pt 1000 stainless steel V4A	1	509-67910-00
Temperature sensor Pt 1000 glass-coated	1	509-67920-00
Pt 1000 clamping system	1	509-63100-00
Software Hei-PC-Control (incl. interface cable)	1	592-50500-00
Accessories for heating bath water and oil 1 l, 2 l, 4 l,		See general catalog

Additional accessories are shown in the general catalog or on our website at www.heidolph-instruments.com.

➤ Technical data

MR Hei-End

Supply power (W)	230 V (50/60 Hz)	or 115 V (50/60 Hz)
Operational voltage	100 - 240 V (50/60 Hz)	
Power consumption (W)	825	625
Protection class (DIN EN 60529)	IP 32	
Drive	EC-motor	
Operating mode	continuous	
Speed range (rpm)	30 - 1.400	
Display	digital	
Speed accuracy (%)	±1	
Analog / digital interface	yes	
Heating power (W)	800	600
Hotplate temperature (° C)	20 - 300	
Medium temperature max. (° C)	300	
Accuracy temperature setting (° C)	±1	
Temperature control	Micro controller	
External sensor	Pt 1000	
Temperature accuracy with external temperature sensor* (° C)	±0,2	
Accuracy hotplate* (° C)	±5	
Safety circuit Hotplate (° C)	10 - 25 over nominal temperature	
Stirring capacity, max (H ₂ O) (l)	20	
Load capacity, max. (kg)	25	
Plate diameter (Ø) (mm)	145	
Weight (kg)	2.6	
Dimensions (l x w x h) (mm)	277 x 173 x 94	

* Accuracy determined with following parameters: 800 ml water in 1 l beaker, form H according to DIN 12 331; temperatur 50 °C; magnetic stirring bars 40 mm; speed 600 rpm; sensor depth 60 mm.

➤ Contact / Technical Service

Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please get in touch with the following address.

For repairs please call Heidolph Instruments or your local authorized Heidolph Instruments distributor.

Warning: Danger of poisoning!

Contaminated units can lead to severe injury or death of our employees!

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact substance
- take adequate protective action towards our parts receiving and service personnel
- mark the pack in accordance with Ordinance on Hazardous Substances

Included at the end of this manual you will find a „Confirmation of condition of unit“.

→ Kindly copy and fill in this form and submit it prior to shipping the unit for repair.

Contact details:

Heidolph Instruments Germany

Heidolph Instruments GmbH & Co. KG
Technical Service
Walpersdorfer Str. 12
D-91126 Schwabach / Germany

Tel.: +49 - 9122 - 9920-74
Fax: +49 - 9122 - 9920-84

E-Mail: service@heidolph.de
www.heidolph.de

Heidolph Instruments North America

Phone: 1-866-650-9604
E-mail: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

Phone: 01799 - 5133-20
E-mail: service@radleys.co.uk
www.heidolph-instruments.co.uk

All other countries

You will find contact details of your local Heidolph distributor at www.heidolph.com

➤ Warranty

Heidolph Instruments provides a 3-year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online (www.heidolph.com). The warranty is valid from the point of registration. The serial number is also valid without registering.

The warranty covers part and manufacturing defects.

In the case of a part or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damages incurred as a result of improper handling or transport.

Warranty claim?

- ➔ Please inform Heidolph Instruments or your local Heidolph distributor should you wish to make a warranty claim.

➤ Confirmation of condition of unit

→ In the case of repair, copy and complete the Confirmation of condition of unit and send it to Heidolph Instruments.

1. Details about the unit

Product number _____

Serial number _____

Reason for repair _____

2. Has the device been cleaned, decontaminated/sterilized?

Yes _____ No _____

3. Is the unit in a condition which does not represent any health threats for the staff of our service department?

Yes _____ No _____

If not, which substances has the unit come into contact with?

4. Legally binding declaration

The customer is aware of being legally liable to Heidolph Instruments for any damages arising from incomplete and incorrect information.

Date _____ Signature _____

Company stamp

Please note

The shipper is responsible for the return of the goods in well packed condition, suitable for the mode of transport.

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This is a translation of the original operating manual. We reserve the right to make any technical modifications without prior notification. This documentation is not subject to revision; the latest version of this documentation can be found on our homepage. The copyright for all texts and images is held by Heidolph Instruments GmbH & Co. KG.
© Heidolph Instruments GmbH & Co. KG