# Intelligent Evaporation

Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

## Operating Manual

Page 2

Rotary Evaporator Hei-VAP Precision  
Operating manual must be read before initial start-up.  
Please follow the safety instructions provided.  
Please keep for future reference.

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<b>Cleaning and Maintenance</b>	<b>42</b>
Cleaning.....	42
Maintenance .....	43
<b>Troubleshooting</b>	<b>44</b>
Malfunctions and resolution.....	44
Overheat protection: reset.....	49
<b>Assembly</b>	<b>50</b>
Electrical connections.....	50
Connexions Electriques.....	51
Mount device .....	52
Assembly parts.....	52
Transportation lock .....	53
Heating bath .....	54
Vapor tube.....	54
Condenser .....	56
Inlet tubes.....	61
Vapor temperature sensor, AUTOaccurate sensor.....	62
Evaporating flask angle and depth .....	62
Control panel .....	64
Connecting peripheral devices .....	65
Connect cooling water .....	65
Mount tubes.....	65
Connect vacuum.....	66
Additional Accessories .....	73
Guard hood.....	73
Safety shield .....	73
<b>Disassembly and Storage</b>	<b>74</b>
Dismantling, transportation and storage .....	74
Dismantling .....	74
Transportation and storage .....	75
<b>Accessories and Spare Parts</b>	<b>76</b>
Scope of delivery.....	76
Accessories.....	76
<b>Attachments</b>	<b>77</b>
Technical data .....	77
Solvent data .....	78
<b>Service</b>	<b>79</b>
Contact / Technical Service .....	79
Warranty .....	80
Confirmation of condition of unit.....	81
Déclaration de non-opposition .....	82

### > Version and variants

#### Version

This manual describes the function, operation and maintenance of the Hei-VAP Precision evaporator.

Version	Alteration date
1.3	01/2018

#### Variants

The units are available in different configurations. Certain features and functions apply to specific configurations only. The varieties are documented in this manual.

### > About this manual

This operating manual is a component of the device described here.

- Please read this manual carefully and obey all safety and warning notices.
- Ensure that every operator reads this manual.
- Ensure that this manual is accessible for every operator.
- Pass on the operating manual to the subsequent owner.

For the current version of this operating manual in pdf format, please go to:

<http://www.heidolph.com/Support/Operation Manuals/Rotary Evaporators/Hei-VAP Series/Hei-VAP Precision>.

In addition, please observe the regional regulations.

## ➤ Terms and signs used

In this manual you will find standardized terms and signs to warn you of possible dangers or give important hints. Please take special notice to these signs to avoid personal injury or damage to the appliance.

The following terms and type of signs (basic drawing elements) are used:

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This device complies with the following EC-Directives\*:

- 2006/42/EC Machinery Directive
- 2014/30/EU Electromagnetic Compatibility Directive

\*Also see attachment „EU Declaration of Conformity“.

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The device has been constructed according to state-of-the-art technology and recognized safety regulations. However, risks may still arise during installation, operation and maintenance./

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The device is suitable for the following use:/  
L'appareil est appropriée pour:

- Distillation to separate substances/  
Distillation, évaporation
- Concentration / Concentration
- Crystallization / Crystallisation
- Powder drying / Séchage de poudre

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Any use which deviates from the device's intended use is considered to be improper. The manufacturer does not accept liability for any damages resulting from non-permitted uses. The risk is carried by the operator alone./

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### Personnel qualification / Qualification du personnel

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- If the equipment is used in corrosive atmospheres, the service life of the equipment will decrease based on concentration, volumes, and frequency of exposures to these corrosive materials, for example concentrated Hydrochloric Acid (HCl).
- If the equipment is used in corrosive atmospheres, the service life of the equipment will decrease based on concentration, volumes, and frequency of exposures to these corrosive materials, for example concentrated Hydrochloric Acid (HCl)./  
Lorsque l'appareil est utilisé dans une atmosphère corrosive, sa durée de vie va diminuer en fonction de la concentration, du volume et de la fréquence d'exposition aux matériaux corrosifs (par exemple de l'acide chlorhydrique HCl concentré).

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## Modifications to the device / Modifications de l'appareil

- No unauthorized changes may be made to the unit.
- No parts may be used which have not been approved by the manufacturer.
- Unauthorized changes result in the EC Declaration of Conformity losing its validity, and the appliance may no longer be operated.
- The manufacturer is not liable for any damage, danger or injuries that result from unauthorized changes or from operating the unit other than as described in this manual.

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## ➤ Connect/disconnect power cord

### Connect power cord

The power cord comes with a three-wire plug and a recess on the bottom side.

Unit plug

Power cord

The socket is located on the back side of the unit.

- ✓ The device is switched OFF.
- Connect plug to the socket with the recess facing down.
- Push plug up against socket until securely attached.
- Connect plug to power socket.

Unit socket

### Disconnect power cord

- Disconnect plug from power socket first.
- Then disconnect plug from socket on unit.

## ➤ Switch unit on and off

The ON/OFF switch is situated on the left side of the device beneath the control panel.

### Switch unit on

- Press the ON/OFF switch on the right side.
- ✓ The display switches on and shows set parameters.
- ✓ The unit is switched on.

ON/OFF switch: ON

### Switch unit off

- Press the ON/OFF switch on the left side.
- ✓ The display switches off.
- ✓ The unit powers off.

## ➤ Heating bath and lift

You can fill the heating bath with various heating bath fluids:

- Tap water (no deionized or distilled water)
- Water soluble polyethylene glycol
- Oil with low viscosity (40 cP) and flashpoint > 285°C.

### **Danger/ Danger Risk of explosion! / Danger d'explosion !**

Heating bath fluid oil with a flashpoint < 285 °C might ignite and lead to uncontrollable thermic reactions and explosion./

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The marks inside heating the bath show the minimum and maximum level for heating bath fluid with evaporating flask completely immersed.

Evaporating flask is immersed by moving up and down via the heating bath lift.

When heating bath lift is moved downwards the evaporating flask displaces already filled in fluid.

To prevent heating bath from overflow:

- Move down the lift with evaporating flask mounted before filling bath.
- Fill in heating bath fluid till the maximum level mark is reached.


Do not let the evaporating flask touch the bottom or the edge of the heating bath.

(To adjust the evaporating flask depth and flask angle please view chapter "Assembly", section "Evaporating flask: Adjust Evaporating flask angle and depth".)

### **Immerse flask: Hand lift**

- ➔ Push down hand lift lever.
- ➔ Move pushed hand lift lever to the left to lower or to the right to raise lift.
- ➔ Release hand lift at the desired position.
- ✓ The hand lift locks into place.

Hand lift lever: upper lift position

### **Immerse flask: Motor lift**

- ✓ Unit is switched on.
- ➔ Push lift key down or up to move motor lift.
- ➔ Release lift key at desired position.
- ✓ The motor lift stops.

Lift keys: up, down

## Start-Up

### Fill heating bath

- ✓ Evaporating flask is immersed.
- Fill in heating bath fluid till maximum level mark is reached.
- ✓ The level of fluid is between the minimum and the maximum mark.

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The heating bath plug is located on the right side in the back of the heating bath.


- ✓ Unit is switched off.
- ✓ Heating bath is locked into place at the base of the unit.

- Connect 7-pole socket to plug of the heating bath with the recess facing up.
- Fix with knurled nut turning it clockwise.

Knurled nut

Plug heating bath

Recess

Feather

### Empty heating bath

---

**Warning / Avertissement:**

**Risk of scalding! Risk of slipping! / Risque de brûlures ! Risque de glissade !**

When emptying the heating bath, you may burn yourself or slip on spilled fluid./

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- 
- ✓ Evaporating flask is removed.
  - Turn knurled nut counter clockwise to unlock it.
  - Remove socket from heating bath.
  - Empty heating bath.
  - Lift bath carefully at both handles and empty.

## ➤ Menu navigation

The basic steps of operating the rotary evaporator Hei-VAP Precision will be explained in this chapter.

The parameters for the distillation are set and monitored via the control panel.

The set parameters are stored. The last used parameters will appear after a restart.

The actual value display is always the start screen. All current actual values are displayed as large numbers and the target values as small numbers.

The programs of the top menu bar are highlighted by turning the select knob Hei-GUIDE. The Hei-GUIDE walks through the programs going from left to right. An „<“ or „>“ shows the direction in which the rotary pushbutton controller can be turned, in order to reach the non-visible menu elements.

The highlighted menu element is activated by pressing the Select Knob Hei-GUIDE and the status screen changes to the respective menu.

The user can change the displayed settings in these menus.

Turning the Select Knob Hei-GUIDE is sufficient to change between individual entries. The selected entry is always highlighted.

The selected entry will be called up when pressing the Select Knob Hei-GUIDE. The desired value can be set when turning the Select Knob Hei-GUIDE.

The lower functions in the menu are activated via the grey softkeys.

<b>1</b>	Pushbutton for the button activation of soft key display	<b>5</b>	Pushbutton rotation start/stop (hard key)
<b>2</b>	Menu line of programs	<b>6</b>	Rotary pushbutton controller (Hei-GUIDE)
<b>3</b>	Pushbutton heating bath on/off (hard key)	<b>7</b>	Pushbutton lift down
<b>4</b>	Soft key display functions	<b>8</b>	Pushbutton lift up

### ➤ **Setting the language**

- ✓ The unit is switched on.
  - 1. Activate the menu line function „Setup“ by using the Select Knob Hei-GUIDE (6).
  - 2. The menu opens by pressing the Select Knob Hei-GUIDE (6).
  - 3. Activate the entry „language“ highlighted in blue.  
The currently set language is located next to the „language“ entry.
  - 4. The marking changes to orange after pressing the Select Knob Hei-GUIDE (6) and is activated for setting the parameter.
  - 5. Select the desired language by using the Select Knob Hei-GUIDE (6).
  - 6. The setting is accepted by pressing it again.  
The text will be displayed in the selected language.
  - 7. You will change to the actual value display by selecting the menu entry „accept (OK)“.
- All changes in the menu view are discarded with soft key „Esc“ (1).

### **Context Help**

The context-sensitive help menu function shows detailed information of the just selected menu item (such as „heating bath“, „rotation“...). This provides additional information to aid operator in use of the Hei-VAP Precision controller.

It can be switched on and off in the setup menu. When the context-sensitive help is shown, it hides all current actual values on the right side.

The context help is switched off when shipped.

→ Activation via the „setup“ menu.

- ✓ The unit is switched on.
  - 1. Activate the menu element „Setup“ by using the Select Knob Hei-GUIDE (6).
  - 2. The menu opens by pressing the Select Knob Hei-GUIDE (6).
  - 3. Activate the entry „context help“ highlighted in blue.
  - 4. The marking changes to orange after pressing the Select Knob Hei-GUIDE (6) and activated for setting the parameter.
  - 5. Select the desired status „on/off“ by turning the Select Knob Hei-GUIDE (6).
  - 6. The selection is confirmed by pressing the Select Knob Hei-GUIDE (6) again. Context help will be displayed.
  - 7. You will change to the actual value display by selecting the menu entry „accept (OK)“.
- All changes in the menu view are discarded with softkey „Esc“ (1).

## ➤ Adjust flask rotation

By increasing flask rotation you can increase your evaporation rate. A shorter evaporation time may prevent thermal damage to your sample.

Depending on flask size and heating bath temperature different rotation speeds may be suitable.

### Speed range and setting

Rotation speed can be set from 10 rpm - 280 rpm.

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At high flask rotation speeds bath fluid may splash or spill over and result in a slippery floor and/or severe scalding./

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➔ Test your selected fluid

- with slowly increasing temperatures
- with different evaporating flasks
- with increasing rotation speeds

➔ Note these results and make them accessible for all individuals.

Depending on the oil and age of oil you are utilizing the heating bath temperature may change the viscosity and the volume of the bath fluid.

Keep records of the exact description/manufacturer information and period of use and repeat the documentation in the determined time interval.

## Operation

### Setting the rotation speed (actual value displays of all modes)

- ✓ The unit is switched on.
- ✓ The rotation is switched off (pushbutton Rotation start/stop).
- ✓ Actual value display

The range can be set from 10 - 280 [rpm] in 1 [rpm] increments. 100 [rpm] are set as default when shipping from the factory.

1. Select the entry „rotation“ by turning the Select Knob Hei-GUIDE (6) and call up the selected entry to be changed by pressing the Select Knob Hei-GUIDE (6) again.
2. Turn the Select Knob Hei-GUIDE (6) to set the desired rotation speed. The rotation speed is readjusted synchronously as an active rotation.
3. Pressing it again confirms the change. The entry is exited after 10 seconds in order to rule out accidental adjustments.
4. Activate the Rotation start/stop (5) pushbutton to start the rotation. The set speed appears on the screen.

### ➤ Setting bath temperature

#### Possible temperature range:

- Water bath application 20°C - 100°C
- Oil bath application 20°C - 210°C

With a pre-set temperature above 75 °C to 80 °C we recommend oil or Polyethylenglycol as heating bath fluid.

Set "Oil" as medium if you work with Polyethylenglycol. Please follow temperature specification provided by manufacturer.

#### **Danger/ Danger Risk of explosion! / Danger d'explosion !**

Potential loss of life/ Danger de mort

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#### **Warning / Avertissement: Risk of scalding! / Risque de brûlures !**

If you operate your unit without a safety hood and at high heating bath temperatures there is a risk of scalding./

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Please wear safety clothes (lab coat, gloves, safety glasses) every time you are in near proximity to the unit during operation./

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### **(Actual value displays of all modes)**

- ✓ The heating bath is filled with heating bath fluid.
  - ✓ The rotary evaporator is operational and switched on at the power switch of the base unit.
1. Switch on the heating bath by pressing the pushbutton heating bath on/off (3).  
The illuminated yellow control ring of pushbutton heating bath on/off (3) shows that the heating bath is activated. An icon [heat waves] in the display shows that the bath is in the heating phase (current flows into the heating coils of the bath).
  2. Select the entry „heating bath” by turning the Select Knob Hei-GUIDE (6) and call up the selected entry to be changed by pressing the Select Knob Hei-GUIDE (6) again.
  3. Pressing it again confirms the change.

## **➤ Calculating the vapor temperature**

### **(with optional vapor temperature sensor only)**

If a vapor temperature sensor is connected during the operation:

→ Switch the unit off and on to activate the program part.

The vapor temperature is displayed on the screen of the control panel when the vapor temperature sensor is connected.

## ➤ Process modes

The Hei-VAP Precision forms a unit with the control panel, the vacuum box and the vacuum valve, in which the vacuum is controlled by opening and closing a valve. The vacuum valve is not required when using a speed-controlled pump.

The ventilation valve in the vacuum box is always open when the vacuum control is switched off.

### **Warning/ Avertissement Imploding glass devices!/ Appareils en verre implosifs !**

Serious injuries from shattering glass./  
Risque de blessure grave en cas de bris de glace.

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### SET pressure mode

The vacuum is controlled and maintained at the set pressure in the SET pressure mode.

#### Menu display

The following values can be changed in the menu view:

Heating bath	Setting of the bath temperature. Setting range is between 20 °C and maximum heating bath temperature (MAX TEMP is adjustable in Setup) with a highest range of 210 °C.
Rotation	Setting of rotation speed of the distillation flask. Setting range is between 10 and 280 rpm.
Vacuum	Setting of the absolute pressure in the rotary evaporating system. Setting range is between 1 and safety pressure (SAFETY PRESSURE is adjustable in Setup) with a high limit of 1200 mbar.
Max. power pump	Vacuum pump (Vario models) throughput % to maximize pump power. Setting range is between 10 and 100 %. Reduction of the throughput leads to higher solvent recovery rates.
Hysteresis	The value of the hysteresis indicates how much the pressure may rise (in mbar) before the vacuum valve is opened. Setting range is between 1 and 50 mbar (only by using valve controlled pumps Rotavac valve control, Rotavac valve tec, or house vacuum source).
Timer	Setting of the timer to shut down the process. Setting range is between 1 and 1440 minutes. For continuous operation choose „off“ for Timer setting. The heating bath at the end of a timed run can be set to turn off (Parameter found in Setup).

Selecting the entry by turning the Select Knob Hei-GUIDE (6)

1. Call up the selected entry to be changed by pressing the Select Knob Hei-GUIDE (6) (entry will be highlighted in orange).
2. Set the value by using the Select Knob Hei-GUIDE (6).
3. Pressing it again confirms the change. Softkey „Esc“ (1) terminates the change of the currently processed entry.
4. The new settings are accepted by the „accept“ entry. It is changed to the actual value display.  
All changes in the menu view are discarded with softkey „Esc“ (1).  
The parameters are stored under a name to be selected by using softkey „Store Favorite“ (1).

#### Actual value display

- ✓ The evaporator is in the SET pressure mode

The setting range of the menu entry „vacuum“ can be set between 1 and the safety pressure in [mbar], however, not higher than 1200 mbar.

Rotation speed and heating bath can be changed as described in chapter “Operation”, section “Adjust flask rotation” and section “Setting bath temperature”.

## Vacuum control

1. Select entry „vacuum“ by turning the Select Knob Hei-GUIDE (6) and calling up the selected entry to be changed by pressing the Select Knob Hei-GUIDE (6) again.
2. Set the absolute pressure by using the Select Knob Hei-GUIDE (6).  
The pressure is readjusted synchronously at an active vacuum control.
3. Pressing it again confirms the change.
4. The entry is exited after 10 seconds in order to rule out accidental adjustments.

## Softkey „Start all“ and „Stop all“

A distillation can be started and stopped automatically by using softkey “Start all” (1) and „Stop all” (1).

## AUTOeasy mode

The required vacuum is automatically calculated based on a pressure measurement in the “AUTOeasy”. The boiling pressure is detected and maintained automatically.

„AUTOeasy” is configured for easy handling and for common heating bath fluid and cooling agent temperature. In case the amount of steam produced is insufficient and the evaporation does not start, increase the heating bath temperature. Decrease the cooling agent temperature or heating bath temperature in case of insufficient reclaim of solvent.

## Menu display

The following values can be changed in the menu view:

Heating bath	Setting of the bath temperature. Setting range is between 20 °C and maximum heating bath temperature (MAX TEMP is adjustable in Setup) with a highest range of 210 °C.
Rotation	Setting of rotation speed of the distillation flask. Setting range is between 10 and 280 rpm.
Timer	Setting of the timer to shut down the process. Setting range is between 1 and 1440 minutes. For continuous operation choose „off” for Timer setting. The heating bath at the end of a timed run can be set to turn off (Parameter found in Setup).
Pump type	Selects the connected pump type. The internal control parameters are adjusted according to the connected pump (only by using rpm-controlled pump Rotavac vario control or Rotavac vario tec.)

Selecting the entry by turning the Select Knob Hei-GUIDE (6):

1. Call up the selected entry to be changed by pressing the Select Knob Hei-GUIDE (6) (entry will be highlighted in orange).
2. Set the value by using the Select Knob Hei-GUIDE (6).
3. Pressing it again confirms the change. Softkey „Esc” (1) terminates the change of the currently processed entry.
4. The new settings are accepted by the „accept” entry. It is changed to the actual value display.
  - All changes in the menu view are discarded with softkey „Esc” (1).
  - The parameters are stored under a name to be selected by using softkey “Store Favorite” (1).

### Actual value display

✓ Prerequisite: The evaporator is in the AUTOeasy mode.

Rotation speed and heating bath can be changed as described in chapter „Operating the control panel“.

### Softkey „Start all“ and „Stop all“

A distillation can be started and stopped automatically by using softkey „Start all“ (1) and „Stop all“ (1), as described in 6.14.

## AUTOaccurate mode

The vacuum in the AUTOaccurate mode is controlled by the vapor temperature sensor AUTOaccurate.

This mode offers the following advantages:

- Only low quantities of solvent are suctioned by this vacuum pump.
- The vapor point is detected automatically and supplemented for mixtures.
- Multiple solvent distillation is possible without any manual adjustments.
- Knowledge of the thermodynamic material data is not required

The system controls the vacuum so that a temperature increase occurs on the AUTOaccurate sensor via the produced vapor quantity. The temperature increase refers to the starting temperature during the deactivated process. It is also necessary that a constant temperature (Heidolph Rotacool or RotaChill recommended) is present in the condenser prior to the start of the process. The Hei-VAP Precision automatically checks the temperature consistency and informs the user, if necessary. The temperature difference between the condenser temperature and the temperature at the AUTOaccurate sensor used for controlling the vacuum will be increased when increasing the distillation rate (dT). This results in increasing the distilled solvent quantity per time.

✓ Vapor temperature sensor AUTOaccurate -(569-00040-00) is installed (see chapter 10.2). The cooling medium has run for 15 min (AUTOaccurate sensor shows a constant temperature).

A distillation with optimal recurring data results at low values for the distillation rate (dT) (1-3).

Low values result in a slow distillation and high values in a quick distillation.

- Select the optimal general conditions.
- AUTOaccurate mode can only be run in G3 and G6 glassware configurations.

## Operation

### Menu display

The following values can be changed in the menu view:

Heating bath	Setting of the bath temperature. Setting range is between 20 °C and maximum heating bath temperature (MAX TEMP is adjustable in Setup) with a highest range of 210 °C.
Rotation	Setting of rotation speed of the distillation flask. Setting range is between 10 and 280 rpm.
Distillation rate (dT)	Setting of the distillation rate. Setting range is between 1 and 10. A high numeric value leads to a higher rate of distillation.
Max. power pump	Vacuum pump (Vario models) throughput % to maximize pump power. Setting range is between 10 and 100 %. Reduction of the throughput leads to higher solvent recovery rates. (only by using rpm-controlled pump Rotavac vario control or Rotavac vario tec.)
Hysteresis	The value of the hysteresis indicates how much the pressure may rise (in mbar) before the vacuum valve is opened. Setting range is between 1 and 50 mbar. (only by using value controlled pump Rotavac valve control, Rotavac valve tec, or house vacuum.)
Multiple solvents	If multiple solvents is set to „no“ then the pressure is not reduced further after the detection of the first boiling point. If multiple solvents is set to „yes“ then the pressure is reduced until the user defined final pressure is reached.
Final pressure	Setting of the absolute pressure where the distillation process is terminated. Setting range is between 1 mbar and the safety pressure (see setup for Safety Pressure).
Timer	Setting of the timer to shut down the process. Setting range is between 1 and 1440 minutes. For continuous operation choose „off“ for Timer setting. The heating bath at the end of a timed run can be set to turn off (Parameter found in Setup).

Selecting the entry by turning the Select Knob Hei-GUIDE (6)

1. Call up the selected entry to be changed by pressing the Select Knob Hei-GUIDE (6) (entry will be highlighted in orange).
2. Set the value by using the Select Knob Hei-GUIDE (6).
3. Pressing it again confirms the change. Softkey „Esc“ (1) terminates the change of the currently processed entry.

Multiple solvent distillation = no:

Must be selected if the first vapor point must be calculated and the pressure should not be reduced further; for example, when roughly separating 2 solvents.

The resulting pressure is accepted as target value for the vacuum. The pressure (after reaching the vapor point) cannot be reduced further in the system when the temperature at the AUTOaccurate sensor drops.

Multiple solvent distillation = yes:

This must be selected if the pressure should be reduced further after detecting the first vapor point; for example, for completely distilling solvent mixtures.

## Operation

The entry „final pressure“ can be only selected with the entry „ Multiple distillation = yes“. If the temperature on the AUTOaccurate sensor drops, the pressure in the system is further reduced until

- a) the control temperature is reached again or
- b) the final pressure is reached.

The system responds similar to „Stop all“ in case b.) (see 6.14). The new settings are accepted by the „accept“ entry. It is changed to the actual value display.

All changes in the menu view are discarded with softkey „Esc“ (1).

The parameters are stored under a name to be selected by using softkey „Store Favorite“ (1).

If the distillation will not be terminated after reaching a defined pressure, enter value 1 as final pressure.

✓ Valve-controlled pump is installed.

The setting range can be adjusted in the menu entry „hysteresis“ (range 1 to 50 mbar).

✓ Speed-controlled pump is installed.

The setting range is adjustable from 10 - 100 [%] in 10% increments in the „speed“ menu entry.

### Softkey „Start all“ and „Stop all“

A distillation can be started and stopped automatically by using softkey „Start all“ (1) and „Stop all“ (1), as described in 6.14.

## ➤ Timer functions

The timer function stops the distillation after a set period.

The setting range is adjustable from 1 - 1440 [min] in 1 minute increments in the „timer“ menu entry. Setting „off“ is comparable to 0.

Selecting the entry by turning the Select Knob Hei-GUIDE (6):

1. Select the entry „timer“ by turning the Select Knob Hei-GUIDE (6) and call up the selected entry to be changed by pressing the Select Knob Hei-GUIDE (6) again.
2. Set the desired timer function for the desired time by using control knob Select Knob Hei-GUIDE (6).
3. Pressing it again confirms the change.
 

The time is running backwards at the start with pressing the “Start Vacuum” soft key and “Start All” soft key. The distillation is interrupted after completing the set time.

  - The system will be ventilated.
  - The rotation is switched off.
  - The lift is raised (only motor lift models).  
If required change lift settings in “Setup”.

Define the shutoff response of the heating bath in the setup.

1. Select the menu element program „Setup“ by turning the Select Knob Hei-GUIDE (6).
2. Select the entry „Bath Power at End“ by turning the Select Knob Hei-GUIDE (6) and call up the selected entry to be changed by pressing the Select Knob Hei-GUIDE (6) again.  
The heating bath will not be switched off after the termination of the timer at the end of the process „on“ and will shut off if „off“ is selected.
3. Pressing it again confirms the change.
4. Select „accept (OK)“ to change the actual value display.

## ➤ Distillation parameters: Favorites

Up to 30 different data records can be stored.

A data record contains all parameters required for the respective process management based on the modes.

### Storing the data record

- ✓ The vacuum control is not activated.
  - ✓ Parameters are set.
  - ✓ You are in the menu view of a mode (SET pressure, AUTOeasy, AUTOaccurate, gradient)
- Press softkey „Store Favorite“.

---

The entry is set to lowercase letters with element „abc“.  
The entry is set to uppercase letters with element „Abc“.  
The last character on the right side is deleted by briefly pressing the „delete“ element.

---

1. Select the letter by using the Select Knob Hei-GUIDE (6) and combine a name by pressing the Select Knob Hei-GUIDE (6).
2. Accept the name with softkey „Store Name“ (1) or discard the name with „Esc“ (1) and interrupt the storage.

You will return to the menu view of the respective mode after saving.

### Call up a data record (Menu favorites)

- ✓ The vacuum control is not activated.
1. Select the menu element „Favorites“ by turning the Select Knob Hei-GUIDE (6).
  2. Select the program name by turning the Select Knob Hei-GUIDE.
  3. Call up the options by pressing Select Knob Hei-GUIDE; the following options will be available:
    - Apply
    - Rename
    - Copy
    - Delete

### Apply data record

The data record is called up and the stored mode is activated. The actual value display is then changed to the stored parameters on the display.

The process can be started with "Startvacuum" or "Start all".

### Rename the data record

The name can be entered similar to the „save“ function. The system checks for duplicates when naming.

### Copy data record

The data record is copied and saved under a new name. An entry of the new program name will be required. A duplication of the issuance of the program name will be checked.

### Delete data record

Deletes the data record permanently.

## ➤ Gradient function

A pressure ramp can be programmed with the gradient function, with which various pressure values can be accessed. The respective pressures are controlled linearly over time. A pressure drop, i.e. by 50 mbar within 10 minutes, therefore results in a slight increase of the pressure/time curve than the same pressure reduction in 5 minutes. The maximum attainable increase of the pressure/time curve is limited by the pump capacity and the leak rate of the system.

### Programming the pressure ramp (menu display)

The user can set the pressure and time in each line with gradient increments (but not the sum time (total), since this is calculated automatically). The selected element is always highlighted. The displayed gradient increments and the special elements on the right of the menu entries are automatically scrolled upward and downward, since there is a total of 20 gradient increments.

The currently selected entry (a pressure or a time) is called up by pressing the Select Knob Hei-GUIDE and the user can set the desired value by turning the Select Knob Hei-GUIDE. Pressing the Select Knob Hei-GUIDE again will confirm the set value.

Entering and confirming value = „---“ at „Vacuum“ or „Duration“ will delete the respective values. The following value pairs slide upward by one increment.

Last increment with a value of > 0 at duration results in the total processing period.

## Operation

The following values can be changed in the menu view:

Heating bath	Setting of the bath temperature. Setting range is between 20 °C and maximum heating bath temperature (MAX TEMP is adjustable in Setup) with a highest range of 210 °C.
Rotation	Setting of rotation speed of the distillation flask. Setting range is between 10 and 280 rpm.
Max. power pump	Setting the maximum suction power of vacuum pump. Setting range from 10 to 100%. A reduction of the suction power results in a higher return rate of the solvent.
Hysteresis	Setting the pressure difference of the target value at which the vacuum is re-adjusted. Setting range 1 – 50mbar (only by using valve-controlled pump Rotavac valve control or Rotavac valve tec, or House vacuum.)
Enter Pressure/ Time Steps	Changes the sub-menu for programming the pressure-time ramps.

1. The data record can be saved by using softkey „Store Favorite“ (1).
2. Softkey „Esc“ (1) discards the changes and returns to the actual value display.
3. The gradients can be entered in sub-menu „Enter pressure/time steps“.
4. Up to 20 gradient steps can be entered, whereby each gradient step consists of pressure and time.

The pressure and time can be set in each line. However, the sum of time cannot be set, since this is automatically calculated.

The starting point is always at a time value of 0 and the atmospheric pressure.

1. Select the pressure/time step to be changed by turning the Select Knob Hei-GUIDE (6).
2. Confirm the selection by pressing the Select Knob Hei-GUIDE (6).
3. The value can be changed by turning the Select Knob Hei-GUIDE (6) to the left or right.
4. The current setting is accepted by pressing the Select Knob Hei-GUIDE (6).  
The data record can be saved by using softkey „Store Favorite“ (1). Softkey „Esc“ (1) discards the changes and returns to the actual value display.  
A graphic overview of the programmed ramp is displayed by using softkey „Graph“ (1).

The gradient diagram presents the time on the x-axis in [min] and the related pressure on the y-axis in [mbar].

The scaling of the x and y axes is illustrated so that all of the entered times or pressures are illustrated.

The entry „enter basic settings“ returns to the time-unrelated basic values of the heating bath, rotation, hysteresis, and the max. power pump.

## Start pressure ramp

The pressure ramp is started with softkey „Start all“ (1) or „Start vacuum“ (1). The pressure/time segments are processed.

Note: The maximum attainable increase of the pressure/time curve is limited by the pump capacity and the leak rate of the system.

The process is terminated when reaching the last entry. The vacuum control is switched off, the system is ventilated. The lift is raised with a motorized lift. The heating bath is switched off depending on the setting in the setup.

## Calling up the pressure ramp

Pressure ramps can be called up via the Favorites menu.

The programmed pressure ramps can be edited when the process is deactivated by changing the respective parameters (see chapter 6.12.1).

## ➤ Softkeys

### Softkey „Graph“

The pressure curve is displayed in the graphic view of the actual value display. The time is illustrated in [min] on the x axis and the associated pressure in [mbar] on the y axis.

Switching over between „graph“ and „values“ is not possible when the vacuum control unit is switched off.

Only the section for which pressure values are available is shown on the pressure axis. The displayed time range on the x axis is set in the setup. If the set time range is exceeded, the time window will be continuously shifted.

### Softkey „Start all“ and „Stop all“

A distillation can be started and stopped automatically by using softkey „Start all“ (1) and „Stop all“ (1).

- ✓ The unit is operational.
- ✓ The immersion depth of the evaporator flask is set (see chapter 4.5.7).
- ✓ The desired settings are completed.

Caution - Glass breakage! The immersion depth of the evaporator flasks must be set correctly to prevent lowering of flask into bot-tom of heat bath on motor lift models.

### “Start all“

Caution - Glass breakage! The immersion depth of the evaporator flasks must be set correctly to prevent lowering of flask into bot-tom of heat bath on motor lift models.

All parameters of the distillation can be changed during the course of the „Start all“ function.

### “Stop all“

The lift moves upward (only for the motor lift unit), rotation, vacuum control and timing are stopped.

## Overview of all softkeys

Softkey	Function
Store Favorite	Store data record for Favorites (program list).
Esc	Cancel entry and back up on step, no storage of changes
Hold	Vacuum valve closed, vacuum is kept; no ventilation
Start all	Vacuum control on, rotation on, lower motor lift
Stop all	Vacuum control off, rotation off, raise motor lift
Start Vacuum	Vacuum control on
Stop Vacuum	Vacuum control off, ventilate system
Graph	Display of diagram of pressure/time curve The actual curve is displayed in the actual value display. The tabular entries are visualized in the menu program gradient.
Max. Vacuum	Vacuum control on, system is evacuated to the minimum pressure value allowable considering pump and air leaks.
Continue	Vacuum control on (after Hold)
Values	The actual value screen or gradient increments screen is displayed.
Accept	The current actual value of the pressure is accepted as a new set value.
List	Display of all lists of the stored Favorites.
Scroll Help	Scrolling within the context help to visualize hidden text.

## ➤ More functions

### Calibrating the vacuum sensor

The vacuum sensor can be calibrated.

- ✓ Vacuum gauge is available.
- ✓ The unit is connected.
- 1. Interconnect the vacuum gauge near the condenser (for example, with a Y-piece).
- 2. Switch on the vacuum pump.
- 3. Change to the setup menu and select calibration p-sensor.
- 4. Confirm with „Yes“.
- 5. Enter the atmospheric pressure calculated by an external pressure gauge and confirm by pressing the Select Knob Hei-GUIDE.
- 6. The vacuum control now automatically reduces the pressure in the system. Please wait until the pressure value is stable.
- 7. Enter the displayed pressure calculated by an external pressure gauge and confirm by pressing the Select Knob Hei-GUIDE.
- ✓ The calibration is completed.

## Connecting to a PC

The vacuum box can be connected to a PC via USB.

1. Connect the vacuum box to a PC using an USB cable.
2. Configure Hyper Terminal on the PC.
3. Interface settings:
  - Bits per second: 115200
  - Data bits: 8
  - Parity: none
  - Stop bits: 1
  - Stream control: none
4. The protocol values are displayed in Hyper Terminal:\*\*\*  
(only available, if rotation is switched on)
  - For example: 20;31.4;25.3;976
  - Rotation = „20“
  - HB-Temperature = „31.4“
  - Ext. Sensor-Temperature = „25.3“
  - Vacuum = „976“

## ➤ Setup

Online Help	Shows Online-Help instead of the actual values in the menu SETpressure, AUTOeasy, AUTOaccurate, Gradient, Favorites and Setup.
Safety pressure	If the system reaches the setting of the Safety Pressure point, the vacuum control and flask rotation is stopped, the evaporator is vented and the motor lift is raised (motor lift only). Setting range is between 900 and 1200 mbar.
Heating Bath Fluid	Allows switching between different heating bath fluid and sets a limit to the maximum heating bath temperature for the fluid for operator safety.
Max. Temp. Heating Bath	Sets the limit of the maximum heating bath temperature which can be set. The chosen bath fluid has an influence on the temperature range. Maximum set range between 20 and 210 °C.
Language	Sets language of the interface.
Bath Power at End	Sets the state of the bath after the process has ended. When setting is „off“ the bath is powered down after the process is terminated.
Brightness Display	Changes brightness of the display.
Calibrate p-sensor	Opens sub-menu for calibration of the pressure sensor.
Error history	Shows the errors and history of occurrence for service data.
Reboot setup	All settings of the setup are put into the state of factory setting when Reboot engaged.
Scale time axis	Sets range of the visible time frame in the graphical „display actual values“. Set range is between 1 and 120 minutes.

## ➤ Menu topology

### Actual value display for all modes

Menu entry	Value range
Vacuum	1 and 1200 mbar; default value: 1000 mbar
Vapor	Vapor temperature measured at the vapor temperature sensor.
Rotation	10 and 290 rpm in 1 rpm increments; default value: 100 rpm
Heating bath	20 °C and the max. temperature of the heating bath (to be assigned in Setup), however, no more than 210 °C
Softkey 1	Start Vacuum
Softkey 2	Start all / Stop all
Softkey 3	Max. Vacuum
Softkey 4	Stop Vacuum

### SET pressure mode

Menu entry	Value range
Heating bath	20 °C and the max. temperature of the heating bath (to be assigned in Setup), however, no more than 210 °C.
Rotation	10 and 290 rpm in 1 rpm increments; default value: 100 rpm
Vacuum	1 and 1200 mbar; default value: 1000 mbar
Hysteresis	1 to 50 mbar; default value: 10 mbar
Timer	1 to 1440 min; in 1 min; default value: 0 min; 0=off
Accept (OK)	
Softkey 1	Store Favorite
Softkey 4	Esc

### AUTOeasy mode

Menu entry	Value range
Heating bath	20 °C and the max. temperature of the heating bath (to be assigned in Setup), however, no more than 210 °C.
Rotation	10 and 280 rpm in 1 rpm increments; default value: 100 rpm
Pump type	"vario control" and "vario tec" ; default value: "variocontrol" Blank out by using vacuum-valve
Timer	1 to 1440 min; in 1 min; default value: 0 min; 0=off
Accept	
Softkey 1	Store Favorite
Softkey 4	Esc

## AUTOaccurate mode

Menu entry	Value range
Heating bath	20 °C and the max. temperature of the heating bath (to be assigned in Setup), however, no more than 210 °C.
Rotation	10 and 280 rpm in 1 rpm increments; default value: 100 rpm
Distillation rate (dT)	<ul style="list-style-type: none"> <li>▪ 1 to 10</li> <li>▪ 1-3 „low“</li> <li>▪ 4-6 „medium“</li> <li>▪ 7-10 „high“</li> </ul>
Hysteresis	Hysteresis value 1 to 50 mbar; default value: 10 mbar
Multiple solvent	yes/no; default value: No
Final pressure	Final pressure 1 mbar ... Safety pressure in 1 mbar; default value: 1000 mbar
Timer	1 to 1440 min; in 1 min; default value: 0 min; 0=off
Accept	
Softkey 1	Store Favorite
Softkey 4	Esc

## Operation

### Gradient

Menu entry	Sub-menu	Element	Value range
Heating bath			20 °C ... max. temp. of the heating bath not higher than 210 °C.
Rotation			10 ... 280 rpm in 1 rpm increments; default value: 100 rpm
Hysteresis			1 to 50 mbar; default value: 10 mbar
Enter pressure/ timesteps	Pressure/Time		
		Step	Program step; 20 program steps consisting of the duration and pressure can be entered.
		Vacuum	1 to 1200 mbar; default value: 1000 mbar
		Duration	Duration of a program step in minutes. The pressure is reduced linear or increased as ramp during this time, depending on which pressure was entered during the previous step.
		Sum of Time (Total time)	Adds the duration of the previous steps. The number specifies the time since the start of the process. The last number equals the total processing time.
	Enter basic Settings	Jump back to main menu „gradient“	
	Accept		
	Store Favorite > store under		
		abc	abc Sets the entry to lowercase letters.
		ABC	ABC Sets the entry to uppercase letters.
		Delete	Briefly pressing the Select Knob Hei-GUIDE deletes the last character on the right side, pressing for an extended period (>1.5 s) deletes all characters of a program name.
		[Character]	a b ... z A B ... Z 0 ... 9 . ; ... `
		enter at least one character	
Accept			
Softkey 1	Store Favorite	Softkey 4	Esc
Softkey 3	Graph/Values		

## Operation

### Favorites

Menu entry	Value range / status
(empty)	
[name of Favorite]	<ul style="list-style-type: none"><li>▪ Apply</li><li>▪ Rename</li><li>▪ Copy</li><li>▪ Delete</li></ul>
Back	
Softkey 4	Esc
Cancel entry and back up on step, no storage of changes	

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## Mount and remove evaporating flask

### Mount evaporating flask

- ➔ Insert evaporating flask with Easy-Clip open over the vapor tube.
- ➔ Tilt Easy-Clip over flange of evaporating flask.
- ✓ Easy-Clip is clicked into place.
- ✓ Evaporating flask is safely assembled.

Easy-Clip open

Easy-Clip closed

### Remove evaporating flask

Removing the evaporating flask:

- Heating must be switched off.
- Rotation must be at „0“.
- Vacuum pump must be switched off.
- Lift with the evaporating flask must be in top position (see chapter "Start-Up, Heating bath and lift").
- Unit must be switched off.

**Danger / Danger:**  
**Risk of scalding! / Risque de brûlures !**

If the evaporating flask has not cooled down it may cause severe scalding. /

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- With 1-2 turns counter clockwise unscrew coupling of the flask clamp Easy-Clip.
- Hold evaporating flask with one hand and open the Easy-Clip.
- Remove evaporating flask carefully.

### Mount and remove receiving flask

#### Mount receiving flask

- Hold receiving flask with open end onto the flange of the condenser.
- Slide clamp over flange of condenser and receiving flask.
- Tighten knurled screw of clamp carefully.
- ✓ The receiving flask is mounted.

#### Remove receiving flask

- Hold receiving flask with one hand.
- Loosen knurled screw of clamp.
- Remove clamp from flange.

Knurled screw

Clamp

### ➤ Dispensing sample

Regulate refill of sample during a process through the inlet tube.

- ✓ Glassware is under vacuum.
- Turn inlet valve counter clockwise.
- ✓ Inlet valve is open and sample is drawn into evaporating flask.
- ✓ The further you open the inlet valve, the faster the sample will be drawn in.
- Turn inlet valve clockwise.
- ✓ Refill of sample will be reduced respectively stopped.

Inlet valve

Ventilation cap

### ➤ Ventilation

The system can be ventilated with the ventilation cap manually.

- To ventilate turn ventilation cap slowly counter clockwise.
- ✓ System is ventilated.
- Turn ventilation cap slowly clockwise.
- ✓ Air inlet is closed gradually.

## ➤ Cleaning

The housing parts and all surfaces of the unit allow for cleaning with a damp cloth and if required a mild soap lotion.

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## **Cleaning instructions for heating bath**

Over a shorter or longer period of time rust and deposition of solids may be observed.

Rust can be removed easily with stainless steel polish.

Lime deposits can be removed with standard cleaning agents based on acetic or citric acid.

Please follow the manufacturer's instructions in both cases.

### Cleaning instructions for coated glassware sets

High temperatures may damage the coating of your glassware sets. Do not exceed a maximum temperature of 50° C when washing coated glassware sets.

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Never continue working with a damaged or a worn-off PTFE seal. Otherwise the unit can be damaged accordingly./

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Seals must be checked on a regular basis and must eventually be replaced./

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The unit is maintenance-free. Any necessary repair must be performed by an authorized Heidolph distributor.

Please contact Heidolph Instruments or your local Heidolph Instruments distributor.

## ➤ Malfunctions and resolution

General errors / message	Cause	Remedy
Unit cannot be switched on	Main plug not connected to the main power supply Fuses defective or blown fuse	<ul style="list-style-type: none"> <li>Connect the main plug with the power supply</li> <li>Replace fuses (see chapter "Assembly, Electrical connection")</li> </ul>
No heating function	Main switch is deactivated	<ul style="list-style-type: none"> <li>Activate main switch</li> </ul>
	Fuses defective or blown fuse	<ul style="list-style-type: none"> <li>Replace fuses (see chapter "Assembly, Electrical connection")</li> </ul>
	Main switch is defective	<ul style="list-style-type: none"> <li>Contact service</li> </ul>
	The connecting cable of the heating bath is not connected	<ul style="list-style-type: none"> <li>Connect the connecting cable of the heating bath</li> </ul>
	The heater of the heating bath is defective	<ul style="list-style-type: none"> <li>Contact service</li> </ul>
	The overheat protection is on	<ul style="list-style-type: none"> <li>If the medium was in the heating bath: Let the heating bath cool and reset overheat protection (see this chapter, section "Overheat protection: reset")</li> <li>If the medium was not in the heating bath: Contact service</li> </ul>
Drive unit is not rotating	Main switch is defective	<ul style="list-style-type: none"> <li>Contact service</li> </ul>
	Speed controller is on the left stop	<ul style="list-style-type: none"> <li>Turn the speed controller to the right</li> </ul>
	Drive unit is defective	<ul style="list-style-type: none"> <li>Contact service</li> </ul>
	Fuses defective or blown fuse	<ul style="list-style-type: none"> <li>Replace fuses (see chapter "Assembly, Electrical connection")</li> </ul>
Motor lift is not functioning	Main switch is deactivated	<ul style="list-style-type: none"> <li>Activate main switch</li> </ul>
	Main switch is defective	<ul style="list-style-type: none"> <li>Contact service</li> </ul>
	Fuses defective or blown fuse	<ul style="list-style-type: none"> <li>Replace fuses (see chapter "Assembly, Electrical connection")</li> </ul>
	Lift is at end stop	<ul style="list-style-type: none"> <li>Activate other arrow key</li> </ul>
	Mechanical system / motor is defective	<ul style="list-style-type: none"> <li>Contact service</li> </ul>
	Height stop is not set accurately	<ul style="list-style-type: none"> <li>Perform height adjustment (see chapter „Assembly, Glassware sets“, section „Install the evaporator flask“)</li> </ul>
No evacuation	Power switch Rotavac vario control or Rotavac vario tec is not activated.	Activate the power switch
	Vacuum valve is defective	
	Fuses defective or blown fuse	Replace fuses (see chapter "Assembly, Electrical connection")

## Troubleshooting

General errors / message	Cause	Remedy
Insufficient vacuum	System leaks	<ul style="list-style-type: none"> <li>Check seals, tubing and connections</li> <li>Check stopcock, grease if necessary</li> </ul>
	Vacuum pump is defective	Observe manufacturer's instructions of the vacuum pump
Unit suddenly switches off	Timer is programmed	Check and switch off the timer, if necessary (see chapter „Operation, Timer functions“)
	Fuses defective or blown fuse	Replace fuses (see chapter „Assembly, Electrical connection“)
	When pressure is detected, the Hei-VAP Precision is emergency stopped.	The Hei-VAP Precision is reset to the starting condition when pressing any key.
	Check the settings in the setup menu.	Calibrate the vacuum sensor (see chapter „Operation“ „More functions“)

Table 7-1: Error codes

Additional conditions Error / message	Reaction
AUTOaccurate sensor is not connected	<ul style="list-style-type: none"> <li>Change to the SET pressure mode</li> <li>The current pressure value is accepted as target value, similar to Hold/Accept in the SET pressure mode</li> <li>AUTOaccurate mode cannot be selected</li> </ul>
Vapor sensor is defective	<ul style="list-style-type: none"> <li>The actual value curve “- - -” is displayed instead of the vapor temperature.</li> </ul>
The vapor sensor is not connected	<ul style="list-style-type: none"> <li>The vapor temperature is not displayed in the actual value display.</li> </ul>
Heater can only be activated at > 0 °C	

Table 7-2: Additional conditions

## Troubleshooting

### Hei-VAP Precision Error / message

The temperature sensor for calculating the vapor temperature is not connected, defective or the temperature is above the measured range.

The vapor temperature sensor is defective or the temperature is below the measured range

The AUTOaccurate sensor is defective or the temperature is below the measured range.

The AUTOaccurate sensor is defective or the temperature is above the measured range.

The transportation protection device is active, please raise the lift and briefly switch the unit off and then on again.

An error occurred in the data bus.

The overheating protection was triggered or the heater is defective

A temperature above the safety temperature was measured in the heating bath. The temperature measurement is defective or the medium filled in was too hot.

### Cause

- Vapor sensor broken
- not connected
- above the temp. range

- Short circuit of the vapor sensor
- below the temp. range

- AUTOaccurate sensor short circuit
- below the temp. range

- AUTOaccurate sensor is broken
- not connected
- above the temp. range

The transportation protection device is still activated.

Cable connection between the control panel and the base unit was disconnected. An error occurred in the data bus.

Heater is defective, i.e. the overheating protection was triggered.

Heating bath safety temperature is exceeded

### Remedy

- Let the vapor temperature sensor cool
- Connect the vapor temperature sensor
- Perform visual check
- Contact service

- Heat the vapor temperature sensor
- Perform visual check
- Contact service

- Heat the AUTOaccurate sensor
- Perform visual check
- Contact service

- Let the AUTOaccurate sensor cool
- Connect the AUTOaccurate sensor
- Perform visual check
- Contact service

- Please raise the lift and briefly switch the unit off and then on again.

- Restore the cable connection between the control panel and the base unit. Switch the power switch off and on again!
- Contact service

- If heating bath fluid was in heating bath: Let the heating bath cool down and reset the maximum temperature limiter (see chapter "Cleaning and maintenance, Maintenance")
- If heating bath fluid was not in the heating bath: Contact service

- Switch off units
- Cool the heating bath medium
- Contact Service if heating bath fluid was not filled in

## Troubleshooting

Hei-VAP Precision Error / message	Cause	Remedy
The heating bath is completely evaporated. Please refill the heating bath medium.	Heating bath is completely evaporated	<ul style="list-style-type: none"> <li>▪ Deactivate the unit</li> <li>▪ Check if the maximum temperature limiter has responded. If yes, reset (see chapter "Cleaning and maintenance, Maintenance")</li> <li>▪ Fill in heat transfer fluid</li> </ul>
An illogical temperature difference was measured between the heating bath sensors. The temperature measurement is defective.	Heating bath sensor variance > 10 K	Contact service
The temperature sensor for the heating bath control is defective. (hardware monitoring unit)	<ul style="list-style-type: none"> <li>▪ Sensor monitoring via hardware</li> <li>▪ Heating bath sensor is broken or short circuited</li> </ul>	Contact service
The temperature sensor for the heating bath control unit is broken or bath medium that is too hot was filled in.	<ul style="list-style-type: none"> <li>▪ Heating bath sensor is broken</li> <li>▪ Heating bath is more than 217 °C</li> </ul>	<ul style="list-style-type: none"> <li>▪ Switch off device</li> <li>▪ Cool down heating bath fluid</li> <li>▪ Contact service, if no heating bath fluid was filled in</li> </ul>
The temperature sensor for the heating bath control unit is defective or bath medium that is too cold was filled in.	<ul style="list-style-type: none"> <li>▪ Heating bath sensor short circuit</li> <li>▪ Heating bath is below 0 °C</li> </ul>	<ul style="list-style-type: none"> <li>▪ Switch off units</li> <li>▪ Cool down heating bath fluid</li> <li>▪ Contact service if heating bath fluid was not filled in</li> </ul>
The lift motor is defective or blocked.	Lift end stop switch is not reached after 10 seconds	Contact service
The set safety pressure was exceeded.	Pressure sensor safety pressure is exceeded	<ul style="list-style-type: none"> <li>▪ Deactivate the unit</li> <li>▪ Check if the pressure sensor has responded. If yes, reset</li> <li>▪ Contact service</li> </ul>
The vacuum box is not connected or defective.	Vacuum box is not connected or defective	<ul style="list-style-type: none"> <li>▪ Check the connection between the vacuum box and the evaporator.</li> <li>▪ Contact service</li> </ul>
The ventilation valve is defective.	Vacuum valve short circuit	<ul style="list-style-type: none"> <li>▪ Check connection between vacuum box and evaporator</li> <li>▪ Contact service</li> </ul>

## Troubleshooting

Hei-VAP Precision Error / message	Cause	Remedy
The vacuum valve is not connected or defective.	<ul style="list-style-type: none"> <li>▪ The vacuum valve is not connected</li> <li>▪ The connection between the vacuum box and the vacuum valve is interrupted</li> </ul>	<ul style="list-style-type: none"> <li>▪ Connect the vacuum valve</li> <li>▪ Check the connection between the vacuum box and the vacuum valve.</li> <li>▪ Contact service</li> </ul>
The vacuum valve is defective.	The vacuum valve is short circuited	Replace the vacuum valve
The vacuum pump is not connected or defective	<ul style="list-style-type: none"> <li>▪ The vacuum pump is not connected</li> <li>▪ The connection between the vacuum box and the vacuum pump is interrupted</li> </ul>	<ul style="list-style-type: none"> <li>▪ Connect the vacuum pump</li> <li>▪ Check the connection between the vacuum box and the vacuum pump</li> </ul>
The rotation drive was overloaded or is defective.	<ul style="list-style-type: none"> <li>▪ Drive unit is defective</li> <li>▪ Motor overload shutdown</li> </ul>	<ul style="list-style-type: none"> <li>▪ Deactivate the unit</li> <li>▪ Stop drive unit for cooling down</li> <li>▪ Contact service</li> </ul>

Table 7-3: Troubleshooting table - Hei-VAP Precision

## Troubleshooting

The rotary evaporator is checked after „power on“ in an initialization phase for connected components. The function scope of the rotary evaporator is based on the connected units.

Missing component	Description and response of the evaporator	Remedy
Vacuum box	<ul style="list-style-type: none"> <li>No vacuum box is connected:</li> <li>Actual value display of the vacuum is not displayed</li> <li>The vacuum entry is not displayed</li> <li>„AUTOeasy“, „AUTOaccurate“ and „Gradient“ modes are not displayed.</li> </ul> <p>Note: An evaporator can be operated with a separate pump stand in this condition.</p>	<ul style="list-style-type: none"> <li>→ Connect component</li> <li>→ Reset via power off/on</li> </ul>
Vacuum valve or vacuum pump	<p>No vacuum valve or vacuum pump is connected.</p> <ul style="list-style-type: none"> <li>Vacuum controller is not active</li> <li>„AUTOeasy“, „AUTOaccurate“ and „Gradient“ modes are not displayed.</li> </ul> <p>Note: An evaporator can be operated with a separate pump stand in this condition.</p>	<ul style="list-style-type: none"> <li>→ Connect component</li> <li>→ Reset via power off/on</li> </ul>
Temperature sensor T boiling	<p>The temperature sensor for calculating the boiling temperature is not connected:</p> <ul style="list-style-type: none"> <li>The vapor temperature is not displayed</li> </ul>	→ Connect component
Temperature sensor AUTOaccurate sensor	<p>The temperature sensor for controlling the AUTOaccurate mode is not connected:</p> <ul style="list-style-type: none"> <li>„AUTOaccurate“ mode is not displayed</li> </ul>	→ Connect component

Table 7-4: Missing Components

## Overheat protection: reset

- Switch unit off. Remove cooled down evaporating flask.
- Disconnect heating bath plug and remove heating bath.
- If necessary empty heating bath and turn it upside down.
- Push the red knob on the bottom of the heating bath with a pointed object.
- ✓ Thermal cut-out is reset.

Reset overheat protection

If you experience a malfunction which can not be resolved, please contact your authorized Heidolph distributor immediately.

## Electrical connections

- Fuses must **only** be replaced by a **professional electrician**.
- **Repairs** of any kind are allowed by **qualified professional electricians only**. Any improper repair can result in a dangerous situation. Contact your local Heidolph Instruments distributor for any repair you may have.

### Warning: Electric shock!



This device must only be connected to a grounded electrical socket.

The unit and voltage must match. The rating plate on the back side of the unit provides all voltage specifics (see picture below).

When shipped, the unit is grounded. In case the original plug is changed, the new plug must have a protective conductor!

If the device is still connected to power while changing the fuse, you might get in contact with live parts.

Prior to replacing the fuse ensure that the unit is switched off and pull the plug from the electric socket. Only use original fuses from Heidolph. After changing a fuse check the device for safe operating condition according to IEC 61010-1.

### Color code for electrical connections:

#### Europe

<b>GREEN/YEL-LOW</b>	PE: Protective conductor ( <b>Earth</b> )
<b>BLUE</b>	N: Neutral conductor
<b>BROWN</b>	P: Phase

#### North America

<b>GREEN</b>	PE: Protective conductor ( <b>Earth</b> )
<b>WHITE</b>	N: Neutral conductor
<b>BLACK</b>	P: Phase

The unit is connected with the power cord supplied to the electrical socket. It plugs into the housing on the rear side of the unit.

For countries utilizing a plug other than the standard one supplied:

- The plug supplied may only be changed by a professional electrician.
- If utilizing an adaptor ensure that it is grounded and approved by the local regulations.
- Wiring, adaptor and plug must at least meet the electrical power stated on the rating plate.

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
## Transportation lock

Transportation lock must be removed before you start. Keep the transportation lock for possible return shipments.

The transportation lock consists of three M5x8 screws and a connection plate with drilled holes to fix the position.

- ✓ Base unit is placed on a stable surface.
- Press and hold lift contact on the lift.

Lift contact

Transportation lock

- Loosen the three screws with allen key supplied and remove plate.
- Release lift contact.
- ✓ Lift drives up automatically.
- Keep all individual parts for subsequent use.


Transportation lock

## Heating bath

### Attach heating bath

- ✓ A movable base plate with an enhanced guide rail is located on the right side of the base unit.
- Place heating bath with gap on the guide rail.
- Move heating bath carefully to the left and right to lock it.

### Move heating bath

It may be necessary to increase the distance between heating bath and drive when bigger evaporating flasks or spacers between the evaporating flask and vapor tube are used.

- ✓ Base unit is placed on a stable surface.
- Allow sufficient space to move heating bath.
- ✓ On the right side the heating bath slides out approximately 20 cm to allow for adjustment.
- Hold heating bath on both handles and slide it to the right until you have reached the desired position.

## Vapor tube



### Assembly parts

- Vapor tube
- Flask clamp Easy-Clip
- Coupling vapor tube
- Clamping sleeve

The flask clamp Easy-Clip is pre-assembled in NS 29. For smaller flasks the Easy-Clip NS 24 is enclosed.

The device is delivered with pre-assembled coupling for the vapor tube and the condenser.


- Hold the locking knob beneath "lock" for a quarter turn.
- Unscrew the coupling of the vapor tube along with the vapor tube on the right side of the drive.
- Unscrew the coupling of the condenser on the left side of the drive.

## Assembly

- Remove tension spring and PTFE seal on the condenser side.

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Especially the sealing lip of the vacuum seal is susceptible to damage in handling. Never operate your system with a damaged or worn-off PTFE seal! Always check the condition of your PTFE seal regularly and if required replace it./

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- ✓ Coupling and Easy-Clip are screwed together.
- Push coupling of the vapor tube along with Easy-Clip on the vapor tube (direction and order as shown under "Assembly parts").
- Push clamping sleeve with the bead to the glass on vapor tube until it audibly locks.
- Insert vapor tube with clamping sleeve first from below to the drive and screw it tight.
- At the very last quarter turn, press locking knob.
- ✓ Coupling locks in and vapor tube is securely mounted.

Vapor tube, clamping sleeve and coupling

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Work with undamaged glass only and check all glassware regularly for any damage./

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Heidolphs coated glassware sets aim to minimize the risk of injury due to glass breakage. However they are only suitable for temperatures up to 80° C.

### **PTFE seal**

- Slide PTFE seal with the labeling "Motor Side" facing drive onto the vapor tube.
- ✓ The PTFE seal is tightened as far as it will go.

PTFE seal

### **Condenser:** **Glassware set G1 and G3**

- Push coupling of glassware over condenser flange.
- Push tension ring over condenser flange.
- Screw condenser to the drive hand-tight.

Tension ring      Coupling

- ✓ Condenser is mounted to the drive.

### Condenser:

#### Glassware set G5 dry ice

- ✓ Condenser is mounted to the drive.
- Screw inlet tube with inlet valve onto opening opposite of vapor tube (see chapter "Assembly, Inlet tubes").
- If using inert gas screw inlet tube with inert gas valve onto opening above inlet valve.

Inert gas valve

- Place red seal ring into notch at top of condenser flange.
- Place white PTFE ring with collar up on to condenser flange.
- Insert the jacket trap into condenser.
- ✓ Due to vacuum inside condenser, the flange is drawn in and the system is tight.
- To protect the dry ice, place black cap on top.

Inlet valve

## Assembly



### Condenser: Glassware set G6 (Condenser for reflux distillation)

Additional assembly parts on the vapor tube:

- PTFE seal ring
- Spacer sleeve
- Additional coupling

- ✓ Vapor tube together with PTFE seal is mounted.
- Insert PTFE seal ring into spacer sleeve.
- Place spacer sleeve including PTFE seal ring with the narrow side facing to vapor tube on condenser side.
- Screw on additional coupling with external thread pointing left.
- Tighten distributor with coupling including tension spring on external thread.

- Insert condenser to the upper joint of distributor.
- Screw inlet tube with inlet valve onto opening opposite vapor tube (see chapter "Assembly, Inlet tubes").


## Assembly

- ✓ Distributor is mounted.
- Turn reflux valve clockwise until O-ring is clearly visible and white plug blocks output to receiving flask.
- ✓ Valve is closed and sample flows back to evaporating flask.
- Turn reflux valve counter clockwise until white plug unblocks the output of receiving flask.
- ✓ Valve is open and sample flows into receiving flask.

(For inlet valve and inlet tube see chapter "Assembly, Inlet tubes".)

The system is not vacuum-tight if valve is opened too far. Open valve only in a way that the white seal ring beneath knurled screw remains visible 0.5 cm wide.

## Assembly

### Condenser support for vertical condenser: Glassware set G3 and G6

#### Assembly parts:

- Support rod
- Boss head with condenser clamp
- Allen key

On the left side of the unit a drilling hole is placed to fix the support rod.

→ Insert support rod with the flattened side downwards into hole.


→ Secure support rod with allen screw using supplied allen key.

→ Slide boss head of cross pole over support rod.

→ Fix cross pole with knurled screw at support rod.

→ Fix the condenser clamp with knurled screw around the condenser.

Knurled screw condenser clamp Boss head

Support rod Condenser clamp Allen key

Hole

Allen screw

Knurled screw boss head

### Condenser support for vertical condenser: Glassware set G5

#### Assembly parts:

- Support rod
- Boss head with tension band
- Allen key

→ Mount condenser support analog to condenser support for G3 and G6.

→ Fix condenser with the tension band.

Tension band Allen key

Knurled screw

Boss head

## Inlet tubes


### Assembly parts:

- 2x Inlet tube
- 1x Inlet valve

On the bottom left side of the condenser a ground opening is located to include stopcock with inlet tube.


### Inlet tube 1

- Shorten supplied PTFE tube so that the end of the tube leads to the evaporating flask and attach it to the stopcock.
- Slide drain disk onto PTFE tube so that it is placed below the vapor opening of condenser.
- ✓ Condensed fluid cannot flow back into evaporating flask through inlet tube.
- Unscrew screw cap from condenser.
- Slide seal ring and screw cap approx. 3 cm onto PTFE tube.
- Insert PTFE tube with seal ahead into inlet opening and tighten screw cap.


### Inlet tube 2

- Unscrew the screw cap including the seal ring of inlet opening.
- Slide seal ring and screw cap approx. 3 cm onto PTFE tube.
- Insert PTFE tube with seal ahead into inlet opening and tighten screw cap.
- ✓ Inlet tubes are mounted.

### Inlet tubes glassware set G1

Openings for inlet tubes are located on the top of the condenser.

## Vapor temperature sensor, AUTOaccurate sensor

**Vapor temperature sensor**  
(Optional accessory)

**AUTOaccurate sensor**  
(Glassware set G3 and G6 only)

Opening for vapor temperature sensor and AUTOaccurate sensor is located at the bottom of the condenser.

- Remove sealing cap of condenser.
- Insert vapor temperature sensor into condenser.
- Tighten vapor temperature sensor with screw cap.
- Plug in connector plug with white mark upwards into receptacle on the left side of unit.

Opening vapor temperature sensor

Connector plug vapor temperature sensor


Screw cap vapor temperature sensor

## Evaporating flask angle and depth

**Evaporating flask angle**

Evaporating flask angle can be adjusted by turning the drive.

- ✓ Rotation is set to 0 rpm.
- ✓ Evaporating flask is mounted (see chapter "Operation", section "Mount evaporating flask").
- Hold condenser with your left hand.
- Slide locking mechanism on lift pillar to the right and hold it against the spring force.
- Move condenser carefully until connected evaporating flask has the correct angle.
- Release locking mechanism.
- If necessary move condenser slightly until it locks in.

Locking

Lift pillar

### Flask depth

A stop bolt is located on the right side of the lift pillar with which you can adjust a constant depth for the evaporating flask.

- ✓ Rotation is set to 0 rpm.
- ✓ Evaporating flask is mounted (see chapter "Operation", section "Mount evaporating flask").
- ✓ Height stop is at top position.

#### **Adjust height without guard hood:**

- Move down lift until it has reached the desired depth (see chapter "Start-Up, Heating bath and lift").
- Turn fixing screw on stop bolt of lift pillar approx. 180° counter clockwise.
- Slide stop bolt via fixing screw downwards as far as it will go.
- Fasten fixing screw clockwise.
- Move the lift upwards.
- ✓ The fixing screw allows to individually adjust the lowest flask position.
- ✓ The lift cannot be lowered beyond this point.

#### **Adjust height with guard hood:**

- Hold guard hood with one hand.
- Move guard hood manually when adjusting the lift.

Fixing screw

Lift contact

Stop bolt

## Control panel

Control panel can be removed from console to be fixed to a wall or fume hood.  
Connection cable is 1 meter long.

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
### **Open cable carrier**

- Grab control panel on the top left side and tip it to the front.
- Avoid pressure to the LCD display.
- ✓ The cable winding is surrounded by a shell.
- Pull shell of cable winding open.
- Wrap cable to the desired length on or off.


### **Position control panel**

- If necessary turn shell 180°.
- Press shell and control panel together until they lock audibly.
- Place control panel on a clean and stable surface or attach it to a wall.

### **Insert control panel into console**

- Turn shell into initial position with the wide side downwards.
- Wrap cable up to the yellow colored cable mark.
- Press shell and control panel together.
- If necessary plug in connecting plug with white colored tag upwards into console.

- Place control panel with an angle of approx. 30° forward into the console.
- Press control panel into console until it audibly locks.

### ➤ Connecting peripheral devices

#### Connect cooling water

##### Assembly parts:

- 2x tube nipples  
(included in scope of delivery)
- 1x tube, inner Ø 7-8mm  
(not included in scope of delivery)


Input and output for cooling water is located on the back of condenser.

- Remove screw cap from input and output of the cooling water.
- Push screw cap over tube nipples so that the black seal is placed inside the cap.
- Push intake and outflow hoses for cooling water through screw cap over one of each tube nipple.
- Tighten screw cap with seal side up to the openings of condenser.
- ✓ Intake and outflow of cooling water are connected.

Tube nipples

Input and output cooling water

#### Mount tubes

To mount tubes a tubing clamp is located on the back of device.

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To connect the vacuum you need:


- a vacuum source
- a vacuum tube ( $\varnothing$  7-8 mm inside)

#### Glassware set G1, G3 and G6

The opening of the vacuum connection is located on the top of condenser.

- Unscrew either one of the screw caps with nipple.
- Push vacuum tubing onto nipple.
- Tighten screw cap with nipple and tubing.
- Connect other end of vacuum tube with vacuum source.
- If necessary connect second vacuum tube to Vacuum sensor

#### Vacuum connection G1, G3 and G6

Screw caps vacuum

#### Glassware set G5

The opening for vacuum connection is located on the top of condenser.

#### Vacuum connection G5


Opening vacuum connection

# Assembly

## Vacuum box and sensor

The support for vacuum box and vacuum sensor is located at the back of device.


Vacuum box

Vacuum sensor

### Connections Vacuum box

#### Connect vacuum box

- Place vacuum box on the guide rails behind the heating bath.
- The vacuum box brackets must be located above the holes.
- Place 2 enclosed screws in the holes with thread through the bracket of the vacuum box.
- Tighten screws.



## Assembly

### Connect Woulff bottle


- Unscrew the threaded fitting from the Woulff bottle.
- Remove the hose clamp ring.
- Place the enclosed PTFE tube section in the cap nut and tighten the cap nut.
- Push the threaded fitting on the PTFE tube.
- Push on the hose clamp ring.
- Connect to the Woulff bottle with the threaded fitting.
- ✓ The vacuum valve and Woulff bottle are connected.
- Attach the connected Woulff bottle and the vacuum valve with the mounting brackets and two screws to the base unit.

## Assembly

### Hei-VAP Precision with Rotavac valve control and vacuum box


<b>1</b>	Condenser	<b>8</b>	Ventilation / inert gas
<b>2</b>	Y-piece	<b>9</b>	Vacuum tube
<b>3</b>	Vacuum valve	<b>10</b>	Vacuum valve cable connection
<b>4</b>	Condensate cooler (optional)	<b>11</b>	USB PC connection
<b>5</b>	Vacuum valve cable	<b>12</b>	Vacuum sensor connection
<b>6</b>	Switch box connection	<b>13</b>	Vacuum sensor
<b>7</b>	Hei-VAP connection	<b>14</b>	Data cable

This installation instruction refers to the current configuration with external vacuum sensor. If you use an older version of this model, please refer to the correct manual.

## Assembly

### Hei-VAP Precision with Rotavac valve control, vacuum box and Wouff bottle


<b>1</b>	Condenser	<b>8</b>	Ventilation / inert gas
<b>2</b>	Wouff bottle	<b>9</b>	Vacuum tube
<b>3</b>	Vacuum valve	<b>10</b>	Vacuum valve cable connection
<b>4</b>	Condensate cooler (optional)	<b>11</b>	USB PC connection
<b>5</b>	Vacuum valve cable	<b>12</b>	Vacuum sensor connection
<b>6</b>	Switch box connection	<b>13</b>	Vacuum sensor
<b>7</b>	Hei-VAP connection	<b>14</b>	Data cable

This installation instruction refers to the current configuration with external vacuum sensor. If you use an older version of this model, please refer to the correct manual.

## Assembly

### Hei-VAP Precision with Rotavac vario tec and vacuum box


<b>1</b>	Condenser	<b>8</b>	Vacuum tube
<b>2</b>	Y-piece	<b>9</b>	Vacuum valve cable connection
<b>3</b>	Condensate cooler (optional)	<b>10</b>	USB PC connection
<b>4</b>	Vacuum valve cable Rotavac vario control / Rotavac vario tec	<b>11</b>	Vacuum sensor connection
<b>5</b>	Switch box connection	<b>12</b>	Vacuum sensor
<b>6</b>	Hei-VAP connection	<b>13</b>	Data cable
<b>7</b>	Ventilation / inert gas		

This installation instruction refers to the current configuration with external vacuum sensor. If you use an older version of this model, please refer to the correct manual.

## Assembly

### Hei-VAP Precision with Rotavac vario tec and vacuum box and Woulff bottle

<b>1</b>	Condenser	<b>8</b>	Vacuum tube
<b>2</b>	Woulff bottle	<b>9</b>	Vacuum valve cable connection
<b>3</b>	Condensate cooler (optional)	<b>10</b>	USB PC connection
<b>4</b>	Vacuum valve cable Rotavac vario control / Rotavac vario tec	<b>11</b>	Vacuum sensor connection
<b>5</b>	Switch box connection	<b>12</b>	Vacuum sensor
<b>6</b>	Hei-VAP connection	<b>13</b>	Data cable
<b>7</b>	Ventilation / inert gas		

This installation instruction refers to the current configuration with external vacuum sensor. If you use an older version of this model, please refer to the correct manual.


## ➤ Additional Accessories


### Guard hood

- Fix carrier guard hood at height stop with two allen screws supplied.
- Release clamping screw and align guard hood while moving it horizontally.
- Tighten clamping screw
- ✓ Guard hood is mounted.
- Open guard hood with handle.

### Safety shield

- Move lift upwards.
- Pull up safety shield to heating bath and place it so that the guard hood can be opened easily, where applicable.
- Turn eccentric screw 180° and fix safety shield to heating bath.

## ➤ Dismantling, transportation and storage

### Dismantling

**Warning:**  
**Risk of scalding!**

Hot glass surfaces and heating bath fluids may result in scalding prior to dismantling.

Ensure that all parts have cooled down to room temperature.

**Risk of poisoning!**

Minor leakages may result in solvents and toxic substances being released into the environment.

Ensure that fluids and solvents are collected in appropriate vessels.

Ensure that the surface of the unit is not contaminated with solvents and hazardous materials.

To clean the unit you may run a distillation with harmless substance.

**Risk of electric shock!**

If the device is still connected to power during dismantling, you might get in contact with live parts.

Prior to dismantling unit ensure that the unit is switched off and the plug is pulled from the electric socket.

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- ✓ The device is switched off.
- Remove all connections to other devices.
- Remove all glass ware and heating bath.
  - Motor lift**
  - Switch device ON.
  - Move lift downwards as far as it will go.
  - Screw transportation lock while holding down lift contact.
  - Switch the device OFF.
- Disconnect plug from power socket and device socket.

## Transportation and storage

- Store and transport the device and its components only if they were emptied, cleaned and dismantled.
- To carry the device wear safety gloves and lift the unit from below.
- Store and transport the device and parts of device in the original packing material or alternatively in an adequate container to prevent damage.
- Seal the package carefully against unauthorized or accidental opening.
- Store the device in a dry and frost-free place.

Improper storage and transportation may result in damages to the system and the mechanical components.

Avoid any kind of shocks during transportation.

## ➤ Scope of delivery

Component	Variant	Quantity	Product number
Hei-VAP Precision glassware set G3*:	Hand lift	1	563-01300-00
Hei-VAP Precision glassware set G3*:	Motor lift	1	564-01300-00
PTFE/FKM seal		1	23-30-01-01-30
Vapor tube		1	514-00000-01
Clamping sleeve		1	23-30-01-05-31
Transportation lock (motor lift)		1	11-300-006-28
Coupling		1	23-09-03-01-03
Tension ring		1	22-03-02-01-05
Flask clamp Easy-Clip NS 29		1	23-30-01-05-29
Flask clamp Easy-Clip NS 24		1	23-30-01-05-57
Joint clamp		1	515-42000-00
Operating Manual English / German		1	01-005-004-23
Warranty registration / Confirmation of condition		1	01-006-002-78
EC Declaration of Conformity		1	01-001-025-06
Power cord		1	depending on country

\* For product numbers in combination with further glassware sets please refer to our general catalogue or to our website at [www.heidolph.com](http://www.heidolph.com).

## ➤ Accessories

Component	Quantity	Product number
Device fuses (115 V and 230 V)	2	11-300-009-40
Vapor temperature sensor	1	569-00030-00
AUTOaccurate-Sensor	1	569-00040-00
Guard hood	1	569-00010-00
Tube: suitable for vacuum and water	1	591-35000-00
PTFE/FFKM seal	1	23-30-01-06-70
Glassware set		see general catalog
Vacuum pumps and chillers type and size vary		see general catalog

Additional accessories are shown in the general catalog or on our website at [www.heidolph.com](http://www.heidolph.com).

## ➤ Technical data

### Hei-VAP Precision

Standard supply voltage	230 V (50/60 Hz)	
Operational voltage (V)	100 - 240	
Supply power (W)	1400	
Protection class (DIN EN 60529)	<ul style="list-style-type: none"> <li>▪ IP 20 device</li> <li>▪ IP 67 cable heating bath</li> </ul>	
Lift distance (mm)	155	
Drive	EC motor	
Rotation speed (rpm)	10 - 280	
Speed control	electronic	
Heating bath	<ul style="list-style-type: none"> <li>▪ inner Ø: 253 mm, outer Ø: 291 mm</li> <li>▪ material: V4A-steel (1.4404) stainless steel X2CrNiMo 17-12-2</li> <li>▪ Volume 4.5 L</li> </ul>	
Heating capacity (W)	1,300	
Temperature range (°C)	20 - 100 H <sub>2</sub> O / 20 - 210 oil	
Temperature control	electronic	
Temperature accuracy Heating bath (°C)	±1	
Safety temperature protection bath	cut off at 5 °C difference to set temperature	
Overheat protection (°C)	250	
Digital display	<ul style="list-style-type: none"> <li>▪ Heating bath temperature (°C)</li> <li>▪ Vapor temperature (°C), (required: vapor temperature sensor)</li> <li>▪ Rotation speed (rpm)</li> </ul>	
Programming ramps	yes	
Condensing surface (cm <sup>2</sup> )	1,400	
Vacuum controller*	integrated	
Timer	yes	
Airborne noise level	< 85 db(A)	
Evaporation rates (L/h) ΔT* = 40°C (continuous operation)	<ul style="list-style-type: none"> <li>▪ Toluene 8.5</li> <li>▪ Acetone 5.8</li> <li>▪ Ethanol 3.5</li> <li>▪ Water 1.2</li> </ul>	
Weight (without glassware) (kg)	17	
<b>Lift</b>	<b>Manual</b>	<b>Motor</b>
Lifting speed (mm/s)	-	25
Dimensions with glassware set G3 (w x d x h) (mm)	739 x 490 x 887	739 x H 420 x 887

\* Only possible in conjunction with vacuum systems

\*\* ΔT = Difference between the heating bath temperature and the vapor temperature

## Solvent data

The graph shows the relationship between the pressure and vapor temperature of a selection of solvents.

The temperature difference between the vapor temperature and the cooling medium should be at 20 K to result in sufficient condensation.

The temperature difference between the heating bath and vapor temperature should be at 20 K to reach a sufficient Distillation rate (dT).

i.e.: Set a vacuum for a vapor point at 40 °C, set the heating bath temperature at 60 °C.


Figure 11-1: Graph

Torr to mbar conversion: [mmHg]  $\approx$  3/4 [mbar]

## ➤ Contact / Technical Service

### Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please contact us.

For repairs please call Heidolph Instruments or your local authorized Heidolph Instruments distributor.

---

### **Warning: / Avertissement**

#### **Danger of poisoning! / Danger de intoxication !**

Contaminated units can lead to severe injury or death of our employees!./

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When shipping items for repair that may have been contaminated by hazardous substances, please:/

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Included at the end of this manual you will find a „Confirmation of condition of unit“.

→ Kindly copy and fill in this form and submit it prior to shipping the unit for repair.

## Contact details:

### Heidolph Instruments Germany

Heidolph Instruments GmbH & Co. KG  
Technical Service  
Walpersdorfer Str. 12  
D-91126 Schwabach / Germany

Tel.: +49 - 9122 - 9920-74  
Fax: +49 - 9122 - 9920-84

E-Mail: [service@heidolph.de](mailto:service@heidolph.de)  
[www.heidolph.de](http://www.heidolph.de)

### Heidolph Instruments North America

Phone: 1-866-650-9604  
E-mail: [service@heidolph.com](mailto:service@heidolph.com)  
[www.heidolphNA.com](http://www.heidolphNA.com)

### Heidolph Instruments United Kingdom

Phone: 01799 - 5133-20  
E-mail: [service@radleys.co.uk](mailto:service@radleys.co.uk)  
[www.heidolph-instruments.co.uk](http://www.heidolph-instruments.co.uk)

### All other countries

You will find contact details of your local Heidolph distributor at [www.heidolph.com](http://www.heidolph.com)

## ➤ Warranty

Heidolph Instruments provides a 3-year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online ([www.heidolph.com](http://www.heidolph.com)). The warranty is valid from the point of registration. The serial number is also valid without registering.

The warranty covers part and manufacturing defects.

In the case of a part or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damages incurred as a result of improper handling or transport.

Warranty claim?

- ➔ Please inform Heidolph Instruments or your local Heidolph distributor should you wish to make a warranty claim.

## ➤ Confirmation of condition of unit

→ In the case of repair, copy and complete the Confirmation of condition of unit and send it to Heidolph Instruments.

### 1. Details about the unit

Product number \_\_\_\_\_

Serial number \_\_\_\_\_

Reason for repair \_\_\_\_\_

\_\_\_\_\_

### 2. Has the device been cleaned, decontaminated/sterilized?

Yes \_\_\_\_\_ No \_\_\_\_\_

### 3. Is the unit in a condition which does not represent any health threats for the staff of our service department?

Yes \_\_\_\_\_ No \_\_\_\_\_

If not, which substances has the unit come into contact with?  
 \_\_\_\_\_  
 \_\_\_\_\_

### 4. Legally binding declaration

The customer is aware of being legally liable to Heidolph Instruments for any damages arising from incomplete and incorrect information.

Date \_\_\_\_\_ Signature \_\_\_\_\_

### Company stamp

#### Please note

The shipper is responsible for the return of the goods in well packed condition, suitable for the mode of transport.

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This is a translation of the original operating manual. We reserve the right to make any technical modifications without prior notification. This documentation is not subject to revision; the latest version of this documentation can be found on our homepage. The copyright for all texts and images is held by Heidolph Instruments GmbH & Co. KG.
© Heidolph Instruments GmbH & Co. KG